




ONE DOSE ONE PATIENT preserved sulcular opening



EXPASYL





"Restorative procedures require imperative precision and perfect preparation, with the respect of marginal periodontium"

Dr. Lasserre France

"My dentist? I recommend him to everyone"

Lana Dr. Lasserre's patient

DIAGNOSIS

TREATMENT PLANNING

> DATA ACQUISITION > PREPARATION & TREATMENT





The right pressure for perfect sulcular opening



Demonstrated performance^(1,2,3) The solution for every clinical situation

EFFICIENT DRYING

and preserved sulcular opening

PRESERVES

Dr Boulic, France

SEALING /

BONDING

the epithelial attachment and gum line

ASTRINGENT EFFECT:

No bleeding

FAST, SIMPLE and cost-affordable protocol

CLINICAL EFFICACY



Dr Gerdolle, Switzerland

Expasyl[™] Exact offers a wide range of applications

CONVENTIONAL & DIGITAL

impression

Class II and V

PRESERVATIVE CARE



IMPLANTOLOGY

⁽²⁾ Bennani V, Inger M, Aarts JM. Comparison of pressure generated by cordless gingival displacement materials. J Prosthet Dent. 2014 Aug;112(2):163-7.

(3) Klettle T, Schulte C. "Force required to displace Pastes and gels for Gingival retraction", 2012 AADR Annual Meeting (Tampa, Florida) -

Poster Session ID 1364.




ExpasyI™ Exact paste provides excellent benefits with little risks. Its viscosity has been carefully calculated to open the sulcus effectively without damaging the epithelial attachment.

Your procedure is both safe and effective. No need to wait for the tissue healing since it has not been damaged.

- > Your patient doesn't suffer

ExpasyI[™] generates pressure that is 1.7 to 9.2 times greater than other retraction paste available on the market(3)



Note: Hemostasyl is included in this study but is not a retraction compound.

Expasyl[™] generates a pressure that is 37 times less than a cord (143 vs 5396 kPa)(1)

The epithelial attachment is preserved (the integrity of the epithelial attachment is no longer guaranteed above 2400 kPa)(2). The risk of gingival recession and bone resorption are reduced

FORCE REQUIRED FOR MOVEMENT WITH RETRACTION GEL AND COMPOUND (IN VITRO)

Your partner for a natural and aesthetic prosthesis

Superior tissue management great impression




ONE DOSE - ONE PATIENT

- The compound opens the sulcus and dries it (blood, saliva and fluids) completely
- Expasyl[™] facilitates access to the preparation limits and emergence profile for perfect impressions



- Avoid cross-contamination
- Cost-effective : no waste of material, up to 3 teeths
- Works with any composite applicator gun



P I alet e-dentisterie

«Like a traditional impression, a digital imprint requires a precise reading of the limits and no bleeding. For this subgingival preparation covering the dark side of the tooth, the gum tissue quickly repositions itself on the margin, making it impossible to read. Expasyl[™] makes it possible to open the sulcus and dry the place for a right impression. Because reading is easier, the procedure is completed more quickly.»

Dr. P. Lalet, E-dentisterie

RETRACTION CORD VS EXPASYL™ EXACT*

- Accessibility to cervical limits without damaging the gingival tissues is essential.
- Retraction cord causes acute injury as laceration and irreversible trauma to the attachment. The risks are recession, loss of epithelial attachment.
- Moreover if the sulcus is damaged, the chemical agents included in the cord penetrate into deep layer cells. This may induce delayed healing, bone loss and peri-implantitis.
- These phenomenom are inadvisable around implant as they reduce the long-term success of the implant project.⁽¹⁾.
- With Expasyl™ Exact, the preparation limits and ermergence profile are fully accessible. By recording these 2 parameters, you obtain a perfectly adapted prosthesis in terms of mechanic and aesthetic expectations.



«[...] the ExpasylTM system, making it possible to deflect the gums without breaking the epithelial attachment and to create an impression for a prosthesis that meets the state of the art.»

Dr. Jacques CHARON, France, Periodontist

* Expasyl™ Exact: Same paste as Expasyl™

(1) New Zealand dentists' use of gingival retraction techniques for fixed prosthodontics and implants" - Azza Al-Ani. Auckland District Health Board, Greenlane Clinical Centre, Auckland, New Zealand. Vincent Bennani, Nicholas P. Chandler, Karl M. Lyons, W. Murray Thomson. School of Dentistry, University of Otago Dunedin, New Zealand - New Zealand Dental Journal, 2010 Sep;106(3):92-96





Fast, simple and cost affordable protocol



3 minutes Protocol* Time saving Money savings

- Placement: easily inserts into the sulcus in 20 seconds. No cannulas required
- Contact time: works in 1-2 minutes max
- Removal: **easily removed** with a gentle air-water spray



Rinse and gently dry the preparation on a healthy periodontium





- **Time savings**: up to 3 minutes per tooth (especially helpful when preparing multiple teeth)



Position the end of the cannula **Slowly extrude**



to create a closed space. the material.





Let Expasyl[™] Exact act for 1-2 min max. The gum tissue whitens, indicating the compression exerted by the material.



Clear the paste with a gentle air water spray with simultaneous surgical aspiration



The sulcus is opened, and the site is completly dry. No injury. No pain. No bleeding. The impression can be taken.

* contact time: 2 min

- ► Expasyl[™] Exact capsule avoids wasted material
- No need for additional anesthesia or hemostasis
- Impression made at the same clinical time

Inject the paste in the sulcus all around the tooth to ensure the best possible impression quality

EXPASYL — EXACT

Expasyl[™] Exact paste* for all clinical cases



Dr Boulic, France

Peripheral preparation of the treated tooth and place of the temporary tooth. Prosthetic limit slightly placed subgingival to respect the emergence profile and the esthetic insertion of the prosthesis in the smile. Temporary tooth placed for 15 days.



Removal of temporary tooth, rinsing tooth preparation and placement of Expasyl[™] Exact to recover emergence profile. Needed an exsanguinated impression for the digital impression technique.



Screenshot of Expasyl[™] Exact during the digital impression phase. The advantage is to let Expasyl[™] Exact act during the pre-op impression (max. 2min). Once Expasyl[™] Exact has been rinsed off, redo the impression exclusively of the tooth concerned. Note: perfect reading of the margins and the deflection of the marginal gingiva

Tooth prepared before applying Expasyl[™] Exact.



Expasyl[™] Exact placed on the concerned tooth. Product placement for 2 minutes to obtain a clear and exact sulcular opening for the rest of the procedure.



After 2 minutes, Expasyl[™] Exact was easily rinsed off with a air/ water syringe, leading to a distinct sulcular opening. Digital impression (with intraoral scanner) of the tooth before placing the crown.

Dr Gerdolle, Switzerland



Dr Baresel, Germany

Putting in place Expasyl[™] Exact on teeth number 11 and 21 to have a clear sulcular opening.

"Insertion of the paste into the sulcus as vertical as possible, with an insertion speed of 2mm/second"



product for 2 minutes to obtain a clean and precise sulcular opening for the rest of the procedure. Digital impression taken during Expasyl[™] Exact's placement.

Expasyl[™] Exact in place

on teeth. Placement of the



After 2 minutes, Expasyl[™] Exact was rinsed off easily with a air/water syringe, leaving a distinct sulcular opening. We can distinguish the "white" part of the gingiva, which represents the dried sulcus.

* Expasy/[™] Exact: Same paste as Expasy/[™]






Photo 1 of the selected zirconia monoblock crown before it was placed. Photo 2 of the crown placed and cemented for aesthetic prothesis.



To finish, placement of a Full Zircone crown with a conventional glass ionomer sealant.





"Almost 10 years ago, we completely digitalized the impression process in our dental practice." Dr. Ingo Baresel is an expert in taking digital impressions and produces a high quality impression with Expasy1™ Exact.



References	Products
EXPASYL™ EXACT (0.3g CAPSULES)	
Ref. 261 011	Expasyl™ Exact x 50
Ref. 261 010	Expasyl™ Exact x 20
Ref. 260 920	Applicator
EXPASYL™ - KIT & APPLICATOR (1g CAPSULES)	
Ref. 294 100	Minikit: 1 Applicator + 6 Capsules + 12 Straight Cannulas (bendable)
Ref. 260 900	Applicator x1
EXPASYL™ - CAPSULES & CANNULAS (1g CAPSULES)	
Ref. 261 030	Capsules x 20
Ref. 294 010	Capsules x 6
Ref. 261 001	Strawberry Capsules x 20
Ref. 261 005	Straight Cannulas (bendable) x 100
Ref. 261 040	Straight Cannulas (bendable) x 40
Ref. 261 015	Pre-bent Cannulas x 100
Ref. 261 045	Pre-bent Cannulas x 40

Not available in all countries. Please contact your dealer.

Class I medical device - CE - For professional dental use only.

Read the instructions in the notice carefully before use. All information that is essential to ensure correct use of these devices is included in the summary of the product characteristics available on the laboratory's website. Manufacturer: Produits Dentaires Pierre Rolland - France. Update: September 2021

PRODUITS DENTAIRES PIERRE ROLLAND SAS | A company of ACTEON Group Zone Industrielle du Phare I 17 av. Gustave Eiffel I 33700 MERIGNAC I FRANCE Tel + 33 (0) 556 34 06 07 | Fax + 33 (0) 556 34 92 92 E-mail : pierre-rolland@acteongroup.com | www.acteongroup.com

