



BDSI 1st Quarter Offers 2012



POLA NIGHT & DAY



SDI Pola day/night hydrogen/carbamide peroxide in 4/10/50 kits

BUY 4 GET 1 FREE
(lowest value item)

POLA PAINT



8% carbamide peroxide gel, in mascara type presentation, simple to use

BUY 1 GET 1 FREE
(limited stock)

RIVA SELF CURE



Buy any 3 boxes of riva self cure/light cure 50 caps **GET 1 FREE**

Buy any 2 boxes of riva luting/luting plus 50 caps **GET 1 FREE**

seT



Dual cured, self etching, self adhering resin cement

Buy 1 box of seT 50 caps **GET 1 FREE**

SDI amalgam specials available please ring

BDSI Instrument Kits


Surgery Instrument Kits



Luxation Instruments

4 Instruments & Stone **£75.00**
7 Instruments & Stone **£125.00**
Individual instruments **£21.00**

Kit 1



3 Tray Setting Contains:

- 3 Aluminium Instrument Trays
- 1 Pack of Tray Liners • 3 Probes
- 3 Tweezers • 3 Dycal Applicators
- 6 Excavators • 1 Box of Mirror Heads
- 3 Mirror Handles • 3 Ball Burnishers
- 1 Spatula • 3 Amalgam Packers
- 3 Flat Plastics • 3 Carvers

Kit 1 - all for £90.00

Kit 2



Forceps Kit Contains:

- 1 Each of BDS items:
- 308 Upper Anterior Roots • 310 Upper Anterior • 314 Upper Premolar • 326 Upper Roots Bayonet • 329 Right Molars • 330 Left Upper Molars • 349 Lower Incisors, Canines & Premolars • 354 Lower Molars • 361 Lower Roots

Kit 2 - all for £99.00

Kit 3



Full Surgery Kit Contains:

Kit 1 + Kit 2 **plus** 1 each of the following items:

- Wide and Narrow Matrix Bands & Holders • Non and Self Aspirating Cartridge Syringes • Flat and Pencil Bur Brushes • Elevators (Kit of 3) • Scissors (curved and straight)

Kit 3 - all for £249.00

Infection Control Offers

New Pico Evo Disinfector

- Fast 35min cycle
- Internal condenser
- Twin spray arms



Offer includes data logger and 10 clip trays

Prestige Ultraclean II washer disinfector dryer



Offer includes printer or data logger and 10 clip trays!

Prestige Century 3 advanced vacuum autoclave



Includes free printer or data logger, full 2 year parts and labour warranty, free instrument kit!

Acuatic S type vacuum autoclave

- 12 and 4 litre models
- Fast cycle times
- Inbuilt data logger
- Long life PTFE Helix device



Please ring for our latest infection control offers!

Handpiece Offers

BDSI EC-20L 1:1 latch type contra. Screw/Knob type also available

Only £55.00
BUY 5 get 1 free



BDSI EC-30BLP 30,000rpm ball bearing PB contra

£99.00
BUY 3 for £275.00



BDSI ES-30A Straight handpiece. 40,000 rpm

£86.90



BIEN AIR BORA LK

Triple pack with light. Fits multiflex couplings

RRP £2035.00 LESS THAN HALF PRICE £975.00

Full Bien Air range now available



BDSI Handpieces Made In Japan

BDSI ME20-M 22,000 rpm motor in 2 or 4 hole

£137.50



BDSI TCP 450 400,000rpm Japanese push button screw on turbine M4 or B2

£150.00 each
BUY 3 for £399.00



NSK PANA-MAX SERIES

Pana Max push button screw on turbines

TRIPLE PACK £549.00



NSK AS2000 AIR SCALER

NSK AS 2000 Air Scaler Borden or Midwest Fitting

NOW ONLY £299.25



Kavo Handpieces offers available please ring

All prices exclude VAT. Offers valid until 31-03-2012. Subject to availability.

Tel:0845 094 2722 Fax:0845 094 2733 E-mail:sales@bdsi-ltd.co.uk Website:www.ukdentalsupplies.com



BDSI 1st Quarter Offers 2012

INSPECTION LAMP

22W fluorescent magnifying lamp with 5" lens, worktop clamp mount

Great Value
£39.95



INSTRUMENT STERILISATION POUCH SEALING MACHINE

- Dual cutting and sealing action
- Thermostatically controlled heater
- Stainless steel covering
- Inset double -blade cutter for easy and smooth cutting
- Maximum cutting width 25cm
- Adjustable width reel holder

25% off
£221.25



PULP TESTER

- Three speed settings: High, Med, Low
- Automatic turn off

Great Value
Only
£106.25



NEW STYLE ANTHOGRYR CROWN AND BRIDGE REMOVER

6961C

- New ergonomic design
- New 1 click system for all attachments
- Fits all E type air and electric motors



Offer Price **£379.00**

ALGINATE MIXER

- Two speed 120rpm and 250rpm
- Mixing free from bubbles and lumps
- Auto alarm after 20 seconds

Offer Price
Only **£129.00**



AMALGAMATOR

ALMTG7

- Digital LED shows speed and time
- Time adjustable 3 – 60 seconds
- Speed adjustable 2800 opm to 4700 opm

Offer Price
Only **£149.00**



APEX LOCATOR ROOT + ACCESSORIES

Works in dry or wet canal, sound alarm indicating apex, adjustable sound volume.

Only
£211.65



ENDO MOTOR

Can store up to 10 programs for speed and torque. Speed range 100-650rpm. Torque range 0.6-5.2Ncm. Auto reverse function when torque exceeds pre set value. Auto stop when file reaches apex.

£335.75



ENDO MOTOR WITH BUILT IN APEX LOCATOR

Brushless electric motor. Same functions as Endo Motor but with Apex Locator. Also includes 1:1 contra angle handpiece.

Only
£505.75



ELITEDENT MS-1

NEW

Unit comes with two self contained medicament feeders and can also be connected to a water line supply. Dual frequency accepts both 25 and 30k inserts. Detachable and autoclavable sleeves. Two position foot control for normal and turbo. Small worktop footprint supplied with 2 inserts.

Only **£399.00**
with additional **FREE** 25k triple bend insert



BONART M1

INSERTS OFFER
BUY 3 GET 1 FREE!



Magnetostrictive Ultrasonic Scaling Unit

Available in either 25KHz or 30KHz

Offer Price of Only **£299.00**
with **FREE** insert

PIEZO SCALER

- Use direct water supply or water bottle
- Autoclavable EMS style handpiece
- Vibration frequency 28kHz
- 7 tips and endo attachment included
- Buy an additional handpiece with the scaler for **half price** **£72.00**

£389.00



All prices exclude VAT. Offers valid until 31-03-2012. Subject to availability.

Tel:0845 094 2722 Fax:0845 094 2733 E-mail:sales@bdsi-ltd.co.uk Website:www.ukdentalsupplies.com