

KAVO Kerr

www.ukdentalsupplies.com/kavo

SMARTtorque Turbine

Kavo Smarttorque S619L / S619C
16 watt power output

Non-Optic
RRP: £542.00

BDSI PRICE: £319.00

Optic
RRP: £677.00

BDSI PRICE: £359.00

BUY 2 GET 1 FREE



EXPERTtorque Turbine

Kavo Experttorque E680L / E680C
18 watt power output

Non-Optic
RRP: £980.00

BDSI PRICE: £465.00

Optic
RRP: £1,047.00

BDSI PRICE: £515.00

QUAD PACK PRICE £1,545.00



MASTERtorque Turbine

One of the quietest high-speed turbines on the market.
Output of up to 23 watt
Direct Stop Technology -
1 second stop.

RRP £1,430.00

BDSI PRICE: £739.00

BUY 4 GET 1 FREE



NSK

www.ukdentalsupplies.com/nsk

Pana Max 2

NSK Pana Max 2
high speed handpiece
for connection directly
on to your existing
tubing.

RRP: £363.00

BDSI PRICE: £325.00

BUY 1 GET 1 FREE



Pana Max 2 PTL/KV

NSK Pana Max 2
high speed handpiece
for connection on to
NSK PTL or Kavo Multiflex
couplings.

RRP: 431.00

BDSI PRICE: £395.00

BUY 1 GET 1 FREE



NEW & IMPROVED MORE POWER

S-Max M900

NSK S Max M900 turbine
non optic Torque /
Standard head.
Available to fit NSK,
Kavo, W&H couplings.
26 watt power output

TRIPLE PACK

NSK OFFER: £1,095.00

BDSI PRICE: £876.00



MK-dent

High Quality Dental High Speed Handpieces Manufactured in Germany

www.ukdentalsupplies.com/mkdent

ECO LINE

MK-dent HE21KL.
Best value triple pack.
24 watt power output
Compatible with NSK,
Kavo, W&H, Sirona
Couplings.
18 Month Warranty

TRIPLE PACK PRICE

£864.45 Optic

£634.95 Non-Optic



CLASSIC LINE

MK-dent HC20KL
Classic Line Power Edition
28 watt power output.
Compatible with Kavo
couplings only.
18 Month Warranty

TRIPLE PACK PRICE

£1,313.25 Optic

£1,093.95 Non-Optic



PRIME LINE

MK-dent HP21L
Vibration free running.
Sound level less than 55dB.
25 watt power output.
Compatible with
NSK, Kavo, W&H,
Sirona Couplings.
24 Month Warranty

TRIPLE PACK PRICE

£1,272.45 Optic



Electric Micromotors - For Greater Power, Flexibility, Control

Upgrade any dental unit from air to electric - no installation required

NSK NLX NANO PORTABLE MICROMOTOR SYSTEM

Equipped with the NXL NANO Brushless Micromotor which is 20% smaller and lighter than its predecessor. 2 year warranty on unit. 1 year on motor.

Forward / Reverse option

Digital Display. Speed adjustment with 2 memory settings.

LED light brightness adjustable.

Speed range: 2,000 - 40,000 rpm

Maximum Torque: 3.4Ncm

NSK RRP: £2,383.00

BDSI PRICE: £1,715.76 PLUS get Nano 95LS FREE



NSK S-Max M95L

1:5 Speed increasing .

Non-Optic

RRP: £774.00

BDSI PRICE: £557.93

Optic

RRP: £980.00

BDSI PRICE: £706.25



NSK Ti-Max X95L

1:5 Speed increasing.

Non-Optic

RRP: £934.00

BDSI PRICE: £673.05

Optic

RRP: £1,186.00

BDSI PRICE: £854.57



NSK M40 XS

Electric Micromotor. Compatible with Bien Air 4VLM & 4VR 400 hoses

Non- Optic

RRP: £644.00

BDSI PRICE: £464.18

Optic

RRP: £1,026.00

BDSI PRICE: £738.72



BUY 3 GET 1 FREE

BUY 3 GET 1 FREE



www.ukdentalsupplies.com/bienair

OPTIMA PORTABLE MICROMOTOR SYSTEM

Equipped with the MCX-LED Brushless Micromotor with 3 year warranty.

2 year warranty on main Optima unit.

Forward / reverse operation for increased flexibility and versatility

SmartLogic control for adjustable speed. Internal Water.

LED light for optimal illumination. Brushless technology for extended service life

Speed range: 1,000-40,000 rpm

Maximum torque: 2.5 Ncm

CALL TODAY FOR A FREE TRIAL!



Bien Air RRP: £1,541.50

Offer Price: £985.00



Micro series handpieces shorter and lighter with 2 year warranty.

CA 1:1L

Dualook 2 optical glass conductors

Speed 600 - 40,000 rpm.

RRP £733.59

BDSI PRICE: £385.00



CA 1:5L

Increasing speed can be used

for caries removal. Dualook

2 optical glass conductors.

Speed 600 - 200,000rpm.

RRP £902.48

BDSI PRICE: £475.00



MC3

Electric Micromotor

Non-Optic

RRP: £736.00

BDSI PRICE: £429.00

Optic

RRP: £920.00

BDSI PRICE: £539.00



Full range of NSK, Bien Air, Kavo and MK Dent available at discounted prices. Call and speak to a member of our team
www.ukdentalsupplies.com | 01226 380400 | sales@bdsi-ltd.co.uk



ME-20M Air Motor 20,000RPM

24 Months Warranty

£150.00



ECN-20L Latch Grip Contra Angle NSK Compatible

12 Months Warranty

BUY 3 @ £49.00



ES-30A Straight 1:1

12 Months Warranty

BUY 3 @ £98.00

KAVO Kerr

Kavo SMARTmatic S20

Contra Angle 1:1
Max speed 40,000rpm.
For CA burs

RRP £208.00

BDSI PRICE: £124.00



Kavo SMARTmatic Endo S81

Contra Angle 8:1
Miniature head for
Endodontic procedures

RRP: £300.00

BDSI PRICE: £179.00



Kavo SMARTmatic S10

Straight Handpiece 1:1
Max speed 40,000rpm
For HP burs

RRP: £200.00

BDSI PRICE: £119.00



NSK

FX-25M

NSK FX-25M 1:1 Direct drive
push button chuck TD.
Longer neck and smaller head.
Max speed 40,000rpm.

RRP: £233.00

BDSI PRICE: £225.00

BUY 1 GET 1 FREE



FX-15M

NSK FX-15M 4:1 reduction
push button chuck TD.
Longer neck and smaller head
Max speed 10,000rpm

RRP: £269.00

BDSI PRICE: £255.00

BUY 1 GET 1 FREE



NAC-EC

NSK NAC-EC Non-optic latch
type contra angle handpiece
for use with CA/RA burs.
Max speed 25,000rpm

RRP £137.00

BDSI PRICE: £79.50

BUY 5 GET 1 FREE



NSK Prophy-Mate Neo



Air powered tooth polishing system.
Includes 6x Flash Pearl Sachets.
To fit NSK, Kavo, Bien Air, W&H, Sirona couplings.

RRP: £994.00

BDSI PRICE: £733.00

Buy
Prophy-Mate
neo and get
FREE second
handpiece
worth £378*

NSK Flash Pearl Powder



100 x 15g Sachets
RRP: £171.30

BDSI PRICE: £123.34

Buy 2 packs
get a Prophy-
Mate neo
nozzle FREE*

Saeyang SP-CE

Saeyang Microtech. Dental Equipment Manufactured in Korea.

Push Button Contra Angle
Lightweight Micro model
weighs just 58 grams
Max speed 40,000rpm

£88.00

BUY 3 @ £79.00



ES-6

Straight 1:1
Easy to open and close collet chuck
Durable Handpiece
Max speed 40,000rpm

£79.00

BUY 2 @ £75.00



Endo A Class LED

Cordless endodontic motor with
LED light. Reciprocating.
9 Memory Settings
Li-Ion Battery
Speed 60rpm - 900rpm

£475.00





Bonart ART E1

Electrosurgery Unit
Includes 7 piece electrode set

£445.00



Bonart ART M3II

25K / 30K Cavitron Compatible Scaler
with direct Water Feed

£349.00 + 1 Insert FREE



Bonart ART M4Pro

25K / 30K Cavitron Compatible
Scaler with Water bottle supply

£495.00 + 2 Inserts FREE



DTE D7

Ultrasonic Piezo Scaler -
Complete with 6 tips and
Satelec compatible handpiece

£240.00

£295.00 (LED model)

DTE



DTE D5

Ultrasonic Piezo Scaler -
Complete with 6 tips and
Satelec compatible handpiece

£164.00

£185.00 (LED model)



Woodpecker UDS-E

Ultrasonic Piezo Scaler -
Complete with 5 tips and
EMS Compatible handpiece

£320.00

£345.00 (LED model)



Woodpecker UDS-L

Ultrasonic Piezo Scaler -
Complete with 6 tips and
EMS compatible handpiece

£229.00

£259.00 (LED model)

Replacement Scaler Handpieces. DTE / Woodpecker. Satelec / EMS Compatible



£89.95 BUY 2 GET 1 FREE



Woodpecker
DTE

G1
GD1

G2
GD2

G3
GD3

G4
GD4

P1
PD1

P3
PD3

Replacement tips

£12.50

BUY 5 GET 1 FREE



Denjoy Pulp Tester

Includes 4 x lip hooks
2 x probes, 1 x lip hook lead

£95.00



Woodpecker DPEX III GOLD

Apex Locator features LCD screen.
97.71% accuracy.

£149.00



NEW Woodpecker X-Cure

1 Second Cure. High Power LED Curing Light.
Output wavelength 385nm-515nm

£295.00



Linea-tacTAC 400M

Digital Amalgamator
Made in Italy

£195.00



All prices exclude VAT. Offers valid until 20/12/2019. Subject to availability, e&o excepted. Offers cannot be combined

BDIA www.ukdentalsupplies.com | 01226 380400 | sales@bdsi-ltd.co.uk

MEMBER

