

Bien Air[®]
Dental

KAVO

Dental Excellence

MK-dent[®]
GERMANY

ND
NAKAMURA DENTAL Tokyo Japan

NSK

ACTEON

ART
BONART CO., LTD.

WOODPECKER

DTE

BDSI

Dental Equipment & Consumables Issue 20



PSPiX²

Call: 01226 380400
ukdentalsupplies.com




BDIA



Contact details:

Telephone 01226 380400
email sales@bdsi-ltd.co.uk
Website ukdentalsupplies.com
Address Unit 14 Claycliffe Business Park
Cannon Way
Barugh Green
Barnsley
S75 1JU

To place an order please call our friendly staff or go online to the above site.

Orders received before 3.00pm will be despatched the same day with FREE DELIVERY*

*Minimum spend £40.00 exc vat



INDEX

3 in 1 syringes & tips	25	Disposable cups	23	NiTi files	23
Acrylic trimmers	22	Electric micromotor	13	NSK	4,5,9,13 14
Acteon	10	Electrosurgery	12	Paper aprons	24
Adhesives	21	Endo files	23	Paper products	24
Air polishers & powder	8	Endo motor	13	Paper points	23
Air scalers	8	Etchants	19	Piezo handpieces	11
Air motors	2-7	Etch tips	19	Piezo scalers	10
Alginates	19	Face masks	24	Piezo tips	11
Alginate mixing machine	14	Flash pearl	8	Pola whitening	26
Amalgam	19	Gigasonic	24	Polishing cups, brushes & paste	22
Amalgam carriers	25	Gigazyme	24	Portable motor	13
Amalgamator	13	Glass ionomers	20	Pulp tester	14
Anthogyr	6, 13	Gloves	24	PXPIX2 scanner and plates	14
Apex locator	13	GP points	23	Riva	21
Autoclaves	14	GP cutter	12	Safety glasses	25
Autoclave pouches & rolls	25	H files	23	Scaler handpieces	11
Barrier covers	25	Handpieces	3-7	Slow speed contras	2-3
Bibs	24	Handpiece oil & nozzles	8	Surgery consumables	21-23
Bien Air	6, 13	High speed turbines	3-7	Stools	27
Blissonic scaler	8	Implant systems	3,5,6 12	Tubing adaptors	25
Bonart	9, 12	Impression materials	19	Teeth whitening	26
Boutique whitening	26	Impression trays	19	Throat pack	24
Burs	22	Instrument trays & racks	15	Tissues	25
Cartridges	4	Instruments	15-18	Titan scalers	8
Cartridge syringes	25	KaVo	3, 13	Towels	24
Cavitron inserts	9	K files	23	Ultrasonic baths	12
Cements & liners	21	Luxators	17	Ultrasonic scalers	9-10
Composites	20	Magnetostriuctive scalers	9	Water distiller	14
Couplings	3-4	Mikrozid wipes	24	Woodpecker piezo scalers	10
Crown & bridge remover	13	Mixing items	23	X ray viewer	14
Curing lights	12	MK Dent	5	X ray sensor covers	25
Desderman	24	Mouthwash	23		
Disinfection products	24	Nakamura	7		



EU Representative



UK Main Dealer



HANDPIECES

KAVO

Dental Excellence

HIGH SPEED SMARTtorque



S609C Standard head 4 hole screw on 4 port spray ceramic bearings rrp £380



S619L Lux Standard head fits Multiflex 460LED 4 port spray ceramic bearings rrp £677 **S619C** Non Lux rrp £490

HIGH SPEED EXPERTtorque



E680L Lux Standard head 23 watts power 4 port spray reduced back suction exchangeable spray microfilter rrp £1047 Also fits Sirona and NSK couplings **E680C** Non Lux rrp £948

HIGH SPEED MASTERtorque



M9000L Lux low height head 30 watts power low noise level. DST - bur stops within a second rrp £1430. Also fits Sirona couplings

COUPLINGS



460E Standard Multiflex non optic coupling rrp £165



460LE Lux multiflex optic coupling rrp £221 **460LED** multiflex coupling rrp £256



465LRN Lux rotary water adjustment with lock settings rrp £424

SLOW SPEED SMARTmatic



S20 1:1 contra no spray 40,000 rpm CA burs 2.35 rrp £208 **S20S** with external spray rrp £226



S10 1:1 straight no spray 40,000 rpm CA burs 2.35 rrp £200 **S10S** with external spray rrp £217

SLOW SPEED EXPERTmatic



E20C 1:1 contra with internal water non lux CA burs 2.35 rrp £661 **E20L** Lux rrp £848



E25C 1:5 speed increasing non lux with internal water FG burs 1.6 rrp 973 **E25L** Lux rrp £1156

SLOW SPEED MASTERmatic



M20L Lux 1:1 contra Lux with internal water. Small head and triple gear system give best vision & quiet operation. Plasmatec coating improves grip 2 year warranty rrp £702

ENDO HANDPIECES



S81 endo handpiece with 8:1 reduction and max 5,000 rpm. For 2.35 dia. files Full rotation rrp 300

SURGICAL SURGmatic



S201L Lux 20:1 reduction external coolant attached to head rrp £1375 **S201XL** with hexagonal clamping



S11C non optic 1:1 40,000 rpm For HP burs 44.5mm rrp £1100 **S11L** Lux Optic 1:1 rrp £1192



S600LED max speed 40,000 rpm max torque 5.5 Ncm rrp £1283

EXPERTsurg unit



Kavo EXPERTsurg LUX Surgical unit RRP £3261

Features

- One touch auto calibration
- Easy to use interface
- Hygienic design
- Guaranteed precision
- High performance 80Ncm torque

Enquire about our latest price and free handpiece offer

Please call for our latest NSK and KaVo prices and offers tel: 01226 380400

HANDPIECES

NSK
HIGH SPEED BASIC



Pana Max2 with reinforced head and more power. Easy cartridge replacement single water spray rrp £363 Screw on. **PAXSUM4** 4 hole **PAXSUB2** 2 hole



Pana Max2 PTL to fit NSK coupling **PAX2-PTL** KV to fit Kavo coupling **PAX2-KV** rrp 431



Pana Max Plus QD 4 hole spray **PAXQDTU** rrp £431

HIGH SPEED ADVANCED



S-Max M Series integrated shaft bearings for reduced vibration.

Access head **SMAX800** optic **SMAX800L** Torque head **SMAX900** optic **SMAX900L**

Also available for Kavo, W&H, non and optic Bien Air and Sirona optic couplings
All non optic S-Max M rrp £678 Optic £728

CARTRIDGES

Pana Max2 canister type **PAX2SU03** rrp £169.40



S-Max & Pana Max Plus rrp £169.40

SX-MU03 Access **SX-SU03** Torque rrp £169.40



Pana Air Ultra Push **NPASU03** **NPATU03** Mach Lite & Mach QD **NMCSU03** **NMCTU03** Phatelus II & CHQD **NCHSU** **NCHTU** rrp £169.40

Pana Air Chuck **NPAS03** **NPAT03** rrp £115.30

COUPLINGS

PTL CL LED rrp £206.50



KCL-LED rrp £325.90



FULL NSK RANGE AVAILABLE

PLEASE RING FOR LATEST OFFERS OR GO ONLINE

SLOW SPEED



NAC-EC 1:1 Latch rrp 137.50



FX-25M 1:1 PB TD rrp £233.50



FX-15M 4:1 redn. PB TD rrp £263.20



FX-65M Straight HP Bur 2.35 dia TD rrp £233.50



S-Max M95L 1:5 increasing max 200,000 rpm Optic rrp £756

SLOW SPEED AIR MOTORS



FX-205M4 25000 rpm motor rrp £365



M205LG M4 LED self generating optics 24000 rpm internal water rrp £959



S-MAX M205 25000 rpm internal water rrp £533

ENDO HANDPIECES



NSK FX-75M mini head with 16:1 reduction. Full rotation rrp £312

IMPLANT & SURGICAL



SGM-ER20i 20:1 Latch ext & int cooling. Miniature head rrp £522
SGMS-ER20i includes depth indicator rrp £779



S-MAX SG20 20:1 PB ext & int cooling 2,000 max rpm stainless steel rrp £951



SGS-ES 1:1 for surgical burs 2.35 dia 40,000 rpm rrp £415



Ti-MAX X SG25L 1:1 direct drive CA burs 40,000 rpm external water PB titanium body with duracoat rrp £865



Ti-Max X DSG20L 20:1 Optic, Double Sealing System, dismantlable titanium duracoated body rrp £1351
DSG20LH with Hexagon chucking system - holds bur firmer rrp £1373
DSG20 non optic rrp £1245



Ti-MAX X SG65L 1:1 Optic HP Burs 40,000 rpm external water titanium body with duracoat rrp £825



Ti-MAX Z-SG45L 45° 1:3 120,000 rpm. Optic FG Burs external water titanium body with duragrip rrp £1079
SG45 non optic rrp £858. **3 yr warranty**

For latest offers please call 01226 380400 or visit www.ukdentalsupplies.com/nsk

HANDPIECES



Surgic Pro



NSK Surgic Pro+ Surgical micromotor implant system

Features

- Compact body & large LCD display
- Advanced irrigation pump
- Memorises 8 implant systems
- Data log & management functions
- 20-2000RPM
- Torque: 5- 80Ncm
- Program adjustment via foot control

A: Non optic set includes S-Max SG20:1 handpiece RRP £3963

B: Optic set includes Ti Max X-SG20L RRP £5309

C: Surgic Pro+ set includes Ti-Max DSG20L and USB port for recording operations rrp £5907

Please ring for latest offer price and free equipment offers

Vario Surg3



NSK VarioSurg3 Piezo ultrasonic surgery system with LED

Features

- Intuitive Piezo feedback
- Large clear LCD screen
- Wide range of tips
- Lightweight slimline handpiece 50% lighter than predecessor
- 20-2000RPM

RRP £7152

Please ring for latest offer price and free equipment offers

Buy both units get extra discount



HIGH SPEED BASIC



HB12 4 hole screw on 18W power 1 hole spray ceramic bearings **£169**



HB21K Fits Kavo Multiflex 4 hole spray 21W power ceramic bearings **£189**
HB21KL with optics **£225**

HIGH SPEED ECO



HE11 2 hole screw on 18W power 1 hole spray ceramic bearings **£179**
HE17 4 hole screw 4 hole spray **£199**



Chrome coating ceramic bearings 24W power. **HE21** available for **Kavo NSK W&H Sirona**. Optic **£339**
 non optic **£249**

HIGH SPEED CLASSIC



HC21K 23W power. Chrome coating, ceramic bearings Kavo connection only
HC21KL optic **£429**
HC21K non optic **£339**

HC20K Power Edition 28W power. Chrome coating, ceramic bearings Kavo connection only
HC20KL optic **£515**
HC20K non optic **£429**

HIGH SPEED PRIME



Titanium coating, ceramic bearings 25Watt power but only 55DB noise vibration free
HP21 available for **Kavo, NSK, W&H Sirona**
 Optic only **£515**

COUPLINGS

Mk-Dent couplings compatible with major brands



QC4012K 2/3 hole for MK Dent and Kavo Multiflex **£75**



QC4014K 4 hole for MK Dent and Kavo Multiflex **£75**



QC4014NT 4 hole for NSK PTL non optic **£75**



QC4014W 4 hole for W&H Roto Quick non optic **£75**



QC4016NT 6 hole LED for NSK PTL LED **£119**



QC5016W 6 hole for W&H Roto Quick optic **£119**



QC6016KW 6 hole for Kavo with optics and water control **£149**



QC6016SW 6 hole for Sirona with optics and water control **£149**

SAVE WHEN BUYING A TURBINE TRIPLE PACK RING FOR DETAILS

HANDPIECES



HIGH SPEED BLACK PEARL



BP-ECO II PB single spray 12W power screw on non optic Borden and Midwest **£156 Triple pack £439**

HIGH SPEED BORA



BORAL 17W power, halogen optic for unifix. **2 year warranty £385**
Triple pack **£1145** inc coupling
BORALED £469 Triple pack **£1395** inc coupling
BORA non optic **£379** Triple pack **£1099** inc coupling

HIGH SPEED TORNADO X



TORNX01 60 year anniversary handpiece - the most reliable handpiece on the market offered with **5 YEAR WARRANTY** LED 18W power. For unifix coupling **£670** Triple pack inc coupling **£1995**

HIGH SPEED COUPLINGS



BAUNI03 Unifix 4 hole non optic coupling **£115**
UNIF-3 3 hole non optic **£115**
BAUNI04L optic coupling **£159**
BAUNI04LW optic coupling with water adj **£199**

SLOW SPEED CONTRA



BACA1 1:1 blue band PB non optic internal water Classic Line **2 year warranty** rrp 624 **our price £335**
BACA1L Optic **our price £390**



BACA15 1:5 red band PB non optic internal water Classic Line **2 year warranty** rrp 810 **our price £435**
BACA15L Optic **our price £479**

ENDO HANDPIECES



CA81 8:1 5,000 rpm max small head size H 9.4mm Dia 8.5mm rrp 648 **our price £465**

MICRO SERIES CONTRA



Micro Series 14% shorter & lighter

BAMCA1L 1:1 blue band PB optic internal water **2 year warranty** rrp £734 **price £390**
BAMCA15L 1:5 red band Optic rrp £902 **our price £480**



BAPM 1:1 straight internal water non optic **2 year warranty** rrp 418 **our price £225**

MICROMOTORS



Bien Air Optima Portable - just connect to air hose and plug into a power source. **No engineer required.** Immediately transform from air to electric and enjoy the quiet more powerful MXC-LED motor 1-40,000 rpm **3 year warranty** rrp 1541
BAOPT01 **£985**



MC3LED with brushes MC3 LK replacement rrp 920 **our price £560** non optic **£449**

MCXLED brushless can be used with Optima rrp 862 **our price £505** non optic **£410**



MX2LED brushless can use with iOptima display rrp 916 **our price £559**



IMPLANT & SURGICAL



CA20 20:1 non optic for implants rrp £1195 **our price £835**
CA20L Micro Series use with MX-LED Micro Series and MXi micromotors rrp £1211 **our price £845**



Bien Air NEW 3rd generation Chiropro console implant unit with **2 year warranty** MXi LED 3rd generation motor with **3 year warranty**. Simplified display, easy to use compact unit. With CA20L 20:1 contra rrp 3420 **£2395** Without 20:1 rrp £2695 **£1895**



Bien Air NEW 3rd generation Chiropro console implant and oral surgery unit with **2 year warranty**. MXi LED Plus 3rd generation motor with **3 year warranty**. Unit alone **£2095**. With CA20L 20:1 contra **£2595**. With CA 1:2.5L **£2695**

HANDPIECES

ND NAKAMURA DENTAL made in TOKYO JAPAN

SLOW SPEED HANDPIECES



EC-20L Latch Contra 20,000 rpm **£65**
ECN-20L Latch contra (12 points) **£65**
EC-30TL 1/4 turn endodontic **£80.55**
EC-30BL ball bearing 30,000 **£79**



EC-50PK Sealed Propphy Snap on **£65**



EC-50PS Sealed Propphy Screw in **£65**



ES-30A Straight 40.000 rpm **£109**
2 rate **£99 each**

SLOW SPEED AIR MOTOR



ME-20M 4 Hole Midwest 20000 rpm
£181.49 2 YEAR WARRANTY

ME-80M 4 Hole Midwest 5000 rpm
£210.60 2 YEAR WARRANTY



ME-20B 2 Hole Borden 20000 rpm
£181.49 2 YEAR WARRANTY

ME-80B 2 Hole Borden 5000 rpm
£210.60 2 YEAR WARRANTY

OIL & NOZZLES



UK manufacture
handpiece cleaning
and lubricating oil
DS01 300ml can
only **£6 6 rate 5.40**
12 rate £5.10



E type **OILN003 £7.05**



Multiflex **HPN01 £5.95**

Quick **HPN02 £5.95**

Pointed **HPN03 £5.95**

HEADS & SHANKS



EH-20L Latch Head 14 points
fits EC-20L **£26.09 3 rate £23.48**

EHN-20L Latch Head 12 points
fits ECN-20L NSK compatible



EH-50PS Screw in **£26.09 3 @ £23.48**



EH-50PK Snap on **£26.09 3 @ £23.48**




RA-10PK Snap on angle for
straight HP **£22.85**




SHS-EC Shank 1:1 **£44.63**

ORTHODONTIC



Vertical reciprocating contra angle for
stripping. For 2.36mm dia attachment
Stroke 1.3mm. Max speed 5000 rpm.
EGN-50VT £175.00



Stripping starter kit -Ring for details

glows™

Economy
Range

GLOWS HANDPIECE WITH SELF GENERATING LED LIGHT



No coupling needed as this fits directly
onto tubing. Economical Glows self
generating LED turbine handpiece.
Standard head push button chuck with
quattro water spray. Easy cartridge
replacement

SMAXPLUSLEDSG 4 hole £229

SMAXPLUSLEDSG B2 2 hole £229



SUPER-MAX LUX M600L Glass fibre
optic rod quattro spray Fits NSK PTL
SMLLUX600L £250. Triple pack £695



SUPER-MAX PLUS SU non optic
quattro spray. Cartridge as above
SMPSU £149. Triple pack £399

QD COUPLINGS



Glows QD-J
Connector **£77.62**
QDJC01 Borden
QDJC02 Midwest

krafit Korea



Small L81mm and light only 58g push
button contra. Very comfortable
budget handpiece
SP-EC £88 3 rate £79



Straight handpiece 40,000 rpm made in
Korea **ES-6 £79**
2 rate **£75 each**

AIR MOTOR KITS Nakamura motor + straight + latch contra - only **£299** 2 yr warranty motor 1 yr handpieces - good value

HANDPIECES

AIR POLISHERS

NSK



NSK Prophylaxis Mate Neo. Handpiece with 60° nozzle head, powder case, replacement nozzle head 80° 5 x15g sachets of Flash Pearl cleaning powder and maintenance kit. Available to fit W&H, Bien Air, Kavo and Sirona couplings **PMN01 £733**



NSK Flash Pearl cleaning powder **FLASHPS** 100x15g sachets **£120.38**
FLASHPB 4x300g bottles **£85.69**



NSK Perio Mate biofilm eraser. 360° rotation. Available 2 42 or 4 hole and also to fit Kavo, W&H, Sirona and Bien Air rrp £1133 our price **£836.14**

MK-dent®
GERMANY



MK Dent Prophylaxis Line Polishers Ergonomically shaped compact and light 360° Rotogrip system. Available in NSK, W&H, Sirona, Bien Air and Kavo Multiflex connections **PR1011 £495**



MK Dent Prophylaxis Powder in mint, strawberry and lime flavours Extra fine grain size 12-14 µm 300 gram bottle **£15**
Save - buy a pack of 4 for **£47.95**



Prophy 600 Air polishing handpiece **PROP600B2** Borden **£195**
PROP600M4 Midwest **£195**
Replacement Tips **£30**

KAVO
Dental Excellence



Kavo Prophylaxis Flex 4 air polisher with 2 grip sleeves 2 cannulas and 2 powder containers **KAVOP4** rrp 1074 our price **£795**

KaVo's PROPHYpearls are spherical shapes and hence provide optimum cleaning power. Flavours: black currant, mint, orange, peach and neutral or mixed. 80 x 15g sachets **£124.50**



ACTEON



Acteon Air N Go Easy adaptable handpiece for both Supra and Sub gingival treatments. Supra nozzle standard 3 others available at extra cost. W&H Kavo Bien Air Sirona and 4 hole midwest fittings. **ANGE01** rrp 1634 our price **£1025**



Sodium bicarbonate based 5 fresh flavours based on 100% natural aromas or essential oils: neutral, cola, raspberry, peppermint and lemon 4 x 250 bottles per pack **ANGP01** rrp 1634 our price **£71.50**



Economical stain and plaque remover Connects to cart tubing 2 or 4 hole versions **PROP300** **£195**

Prophy 600 polishing powder 130g bottle **PLP01** **£9.95**
6 rate £8.95



AIR SCALERS



KaVo Sonicflex 2000N 6000Hz with 1 universal tip **£379**

Kavo No.5 Universal



Kavo No.6 Sickel



Kavo No.7 Perio



Kavo No.8 Perio XL



KaVo genuine tips rrp 98 price **£85**

GK1 Universal



GK2 Interproximal



GK3 Perio



GK4 Longer Perio



Kavo compatible tips **£35**



Titan Blis-sonic Scaler 5 tips Universal, Sickel, Perio, Perio Right, Left & 4 grips **TITANSW** inc coupling **£645**
TITANK screw on **£559**
Individual tips **£60**



NSK AS 2000 in Midwest **AS2000M** or Borden **AS2000B**. Includes 3 tips and tip changer **£335.02**

NSK



NSK S970 Titanium body and 3 level power ring 6200-6400Hz 3 tips + changer **S970** Non Optic **£766.73**
S970L Optic **£895.90**

S1 Universal **AS2000TIP1**



S2 Sickel **AS2000TIP2**



S3 Perio **AS2000TIP3**



All **£55.20** each
Check website or call for NSK offers

glows™ BUDGET SCALER



Economy Scaler with 3 tips and changer **£179** Tips **£25** each
AIRSCALE01 Borden
AIRSCALE02 Mid West

ULTRASONIC SCALERS

ART. BonART



Bonart M3 Magnetostrictive Ultrasonic Scaling Unit in 25K or 30K. LCD display, 3 power modes. Robust, reliable, popular unit - our best seller!

ART325K (25K) ART330K (30K) £349

All Inserts **£46 each BUY 5 Get 1 FREE**



NEW Bonart M4 Magnetostrictive Scaling unit in 25K or 30K. **Portable** with integrated water bottle, 3 operational modes on clear easy read LCD screen.

Supplied with 2 universal inserts

ARTM4PRO25K (25k) £499

ARTM4PRO30K (30k)



Bonart P3 Piezoelectric Scaling Unit LCD display, 3 power modes. Auto gain control.

Supplied with 4 tips, 2 x Universal, 1 x Flat 1 x Slim ARTP3 £395

Quality BonART tips in Swedish steel

All at £30 each buy 5 get 1 free



Bonart Universal insert

ART25INS01 (25k) ART30INS01 (30k)



Bonart Flat insert

ART25INS02 (25k) ART30INS02 (30k)



Bonart Slim insert

ART25INS03 (25k) ART30INS03 (30k)



Bonart Metal Universal insert

ART25INS05 (25k) ART30INS05 (30k)



Bonart Metal Angled Slim insert

ART25INS05 (25k) ART30INS05 (30k)



Deldent Triple Bend insert

CAV25KTB (25k) CAV30KTB (30k)

Quick Disconnects



Female **1/4ADAPT2**
Male **1/4ADAPT3**

£10 each
£18 pair



Bonart MB3 build in Magnetostrictive Scaler with PC board 2" x 3" LCD display, air control and **2 inserts ARTMB3 £425**



Bonart PB3 build in piezo scaler. Set includes 2" x 3" LCD Screen, PC board, air control, Piezo handpiece, 4 tip changers, **4 tips** (Satelec style) and sterilization box. **ARTPB3 £425**



BonART P6 compact Piezoelectric Scaling Unit with Auto gain control to stabilize and adjust power levels to give maximum performance. **Supplied with 4 tips, 2 x Universal, 1 x Flat 1 x Slim ARTP6 £295**



Universal
ARTTIPS1 (SAT)
ARTTIPE1 (EMS)



Flat Chisel
ARTTIPS2 (SAT)
ARTTIPE2 (EMS)



Rounded Flat
ARTTIPS3 (SAT)
ARTTIPE3 (EMS)



10X Interproximal
ARTTIPS10X (SAT)
ARTTIPE10X (EMS)



10P 45 degree
ARTTIPS10P (SAT)
ARTTIPE10P (EMS)



1S Slim Perio
ARTTIPS1S (SAT)
ARTTIPE1S (EMS)

NSK



NSK Varios 970 non optic / Varios 970 Lux versions available. 3 Scaling tips (G4, G6, G82 x Irrigation bottles, 3 x Tip wrench with torque limiter. AC power cord, sterilization case **VAR970** non optic rrp £2205 **£1627.20** **VAR970LED** rrp £2563 **£2269.73**

ukdentalsupplies.com

tel: 01226 380400

ULTRASONIC SCALERS

NSK



NSK Varios Combi Pro – Two oral care systems in one. Combines ultrasonic treatment and powder therapy into a single professional treatment unit that streamlines procedures and saves time and space. Prophy-Handpiece with 60° nozzle, flash pearl (5pcs. × 15g sachets) VA2LUX-Handpiece, 3 Scaling tips and Irrigation bottle (400mL)
VARIOSCP01 £3247.92.

Satelec **ACTEON**



Satelec Newtron P5 XS B LED F62103 B LED technology easy plaque removal Tip colour coding to power output, Bottle water supply **5 Yr unit warranty**
SATP5XSLED rrp £1984 our price **£1535**



Satelec Newtron P5 B LED F61103 B LED technology easy plaque removal Tip colour coding to power output, **5 Yr unit warranty** **SATP5NEWLED**
rrp £1350 our price **£1045**

ACTEON



Satelec Newtron Booster F60204 Uses Newtron handpiece (F12281) non LED. The tip colour coding system associates each tip with each of the 4 power ranges. Set the colour of the power control dial to match the colour of your tip.
5 Year unit warranty
SATNB rrp £1050 our price **£815**

WOODPECKER

Satelec Compatible



Woodpecker DTE-D5 scaler with Satelec style handpiece includes 6 tips **WOOD31 £249**. Also available with LED handpiece **WOOD66 £299**



DTE D7 scaler with Satelec style handpiece includes 5 tips and endo attachment 2 water bottles no mains connection.
WOOD18 £279 Also available with LED handpiece **WOOD66 £329**

WOODPECKER

EMS Compatible



Woodpecker UDS E scaler with EMS style handpiece includes 5 tips and endo attachment 2 water bottles no mains connection. **WOOD17 £325**
Also available with LED handpiece **WOOD54 £375**



Woodpecker UDS-L scaler with EMS style handpiece includes 5 tips and endo attachment **WOOD01 £229.00**

WOODPECKER



NEW Woodpecker PT Master 5 periodontal treatment, scaler and polishing unit with titanium tips which have a very small vibration amplitude and work in a circular motion so your patient feels less discomfort. Complete with 2 LED handpieces, 2 sterilization boxes, 13 tips and 2 endo attachments. Can be plumbed into cart water supply in addition to the bottle.
WOODPT5 £1099

ULTRASONIC SCALERS

Scaler handpieces



Satelec Newtron B LED slim F12609
rrp £500 **our price £395**



Satelec Newtron B LED slim F12905
white F12900 blue rrp £500 **price £395**



Satelec Newtron B slim F12901
rrp £420 **our price £330**



Satelec Newtron F12281
rrp £400 **our price 315**



Satelec Newtron F12200
rrp £400 **our price 210**



Woodpecker DTE Satelec
compatible handpiece
WOOD19 £89.95
Offers available



Woodpecker EMS compatible
handpiece **WOOD02 £89.95**
Offers available



Woodpecker DTE Satelec compatible
Optic handpiece **WOOD60 £119.50**



Woodpecker EMS compatible Optic
handpiece **WOOD59 £119.50**



MK-dent compatible handpieces have
a steel funnel for extra life and include
a universal DLC 40% longer life tip



MK dent Satelec compatible
handpiece non optic
SC21S £135



MK dent EMS compatible
handpiece non optic
SC21E £135



MK dent Satelec compatible
LED 6 ring optic handpiece
SC21SL £199



MK dent EMS compatible
LED 6 ring optic handpiece
SC21EL £199



NSK VA2-LUX-HP Optic rrp £470

our price £416.18

NSK VA2-HP non optic rrp £262.50

our price £193.68

NSK scaler tips available please call
or go online - **check for bogof offer**



Satelec tips available

Styles 1, 2, 3, 1S, 10P, 10X, 10Z
all at **£59**
EX1.EX2,EX3 **£69**



6 Tip Holder
WOOD08
£14.95



Torque Wrench
WOOD07 £15

Woodpecker tips

£12.50 each BUY 5 get 1 free



G1/GD1 To remove calculus on
supragingival neck of teeth

WOOD14 EMS

WOOD22 Satelec & NSK



G2/GD2 to remove supragingival
heavy calculus

WOOD15 EMS

WOOD23 Satelec & NSK



G3/GD3 to remove calculus of
supragingival subgingival
and interdental

WOOD05 EMS

WOOD24 Satelec & NSK



G4/GD4 to remove all supragingival
calculus

WOOD13 EMS

WOOD25 Satelec & NSK



P1/PD1 to remove subgingival calculus

WOOD16 EMS

WOOD20 Satelec & NSK



P3/PD3 to remove subgingival and
deep perio calculus

WOOD04 EMS

WOOD21 Satelec & NSK



SB1/SBD1 Cavity preparation ball
diamond tip for preparation of the
occlusal surface and neck of tooth

WOOD33 EMS

WOOD34 Satelec & NSK



E1/ED1 0.8mm
U files holder
WOOD32 (Sat)
WOOD09 (EMS)
£25

Woodpecker tips buy 5 get 1 free

tel: 01226 380400

SMALL EQUIPMENT

CURING LIGHTS

MK dent DLC black tips

(diamond coated like carbon)

Up to 40% extra life **£49.50**



Tip 1 scaling

ST01ED EMS compatible

ST01SD Satelec compatible



Tip 2 scaling

ST02ED EMS compatible

ST02SD Satelec compatible



Tip 3 scaling

ST03ED EMS compatible

ST03SD Satelec compatible



Tip 6 Periodontal

ST06ED EMS compatible

ST06SD Satelec compatible



BONART OP1 PIEZO SURGICAL UNIT

ARTOP1 £2250

Features

- Powerful precision cutting
- Touch control display
- 6 tip set
- Controlled torque wrench
- Sterilization box

Applications

- Bone surgery
- Sinus lifts
- Periodontal surgery
- Endodontic surgery
- Orthodontics

LED-C

850-1,000mW, spare battery constantly charging. Three modes, 4 pre set times. Full charge 200 x 10 secs.

LED-C £125



X-CURE

Wide spectrum 1 second curing. Operate in low or high power. Latest OLED display. Easy caries detection. 1,000-3,000mW Wavelength 385nm - 515nm **X-CURE £295**



Ultralite 1000E

LED gun type suitable for small hands. High power 1,600mW, built in light checker, 5 secs to cure composite. Full charge 500 x 5 secs.

ULTRA1000E £325



Ultralite 1800E

LED. High power 2,000mWs. Wavelength 390nm - 480nm. Spare battery constantly charging, built in light checker, 3 seconds fast curing composite. Full charge 500 x 3 second.s **ULTRA1800E £395**

RADII PLUS

High intensity 1500mWm, 5 year warranty. One charge gives 1200 10 second cures. Built in radiometer Operate in low or high power.

RADIIP01 £845

BUY ANY 10 BOXES OF SDI RIVA AND GET A RADII PLUS FREE

Radii attachment heads

5600092 Full Arch Bleaching LED kit **£365**

5600259 Orthodontic (Ortho) LED **£179**

5600260 Right angle standard LED **£174**

5600276 Contra angle standard LED **£174**



RADII-CAL

High intensity 1200mWm, 2 year warranty. One charge gives 1200 10 second cures. **RADIIC01 £649**

ELECTRO SURGERY



Bonart E1 3modes Cut, Coagulation and Cut-coagulation, 10 power levels.

Includes 7 tip electrode set

ARTE1 £495

Electrode Set

7 tips

ARTE103 £82.50



SMALL EQUIPMENT

AMALGAMATOR



TAC-400M Italian manufacture. Ten memory settings, speed 4200 rpm. Accepts amalgam and composite capsules. **ALMT01 £195**



SDI ULTRAMAT 2 Mixing frequency of 4600 oscillations per minute 0-16 second microprocesso, quiet low vibration, 5 year warranty, suitable for all types of encapsulated materials. **ALMT03 £445**



SDI ULTRAMAT S Quiet low vibration, 2 year warranty, suitable for all types of encapsulated materials. **ALMTS01 £375**

BUY ANY 7 BOXES OF SDI RIVA AND GET AN SDI ULTRAMAT S FREE

GP CUTTER



Gutta Percha Cutter with 3 tips. Heats to 500°C in 1 second **MOTGPCUT £90**

PORTABLE MOTOR KIT



Set with E type motor. Rechargeable Li-Ion battery lasts 4-8 hours. Motor speed up to 30,000 rpm 12V 25Watt. **K38MSET £249**

ULTRASONIC BATHS



U95 1.95L working capacity. Internal dims. 200 x 150 x 75 basket 180 x 135 x 45 **U95 £295**

U300 2.5L working capacity. Internal dims. 240 x 140 x 100 basket 220 x 125 x 40 **U95 £319**



Bonart Scooba 4.8L capacity. Heater range 10-60°C, timer 1-99 mins, basket and lid **ARTUC1 £495**

NEW



Hygea d3ntal. Compliant to HTM01-05 and ISO standards. 8 Ltr capacity, heated digitally controlled bath with traceability 6 minute rapid cycle. Timer 1-99 mins **HYGD3 £1675**

APEX LOCATOR



Dpex III Gold The clear LCD screen shows the trajectory of the files using different colours. Multi-frequency apical positioning. Rechargeable lithium battery. Works in wet canals. Adjustable viewing angle **WOOD63 £199**

ENDO MOTORS




Endomate TC2 NSK Cordless endodontic set with torque control and auto reverse. 5 program settings. **ENDO-MATETC2 £1115.86**



Endo A Class Cordless LED endodontic set. Torque level real time display. Reciprocating system. Auto stop and reverse feature. Push button miniature head with 9 program memory. Torque 0.5-3.0Ncm, Speed range 200-600rpm. Li-Ion battery life 90 mins/charge **ENDOA01 £475**

AUTOMATIC CROWN & BRIDGE REMOVER



ANTHOGRYR Automatic C&B remover fits E type motor. Set includes 4 hooks and 2 wire tips **ANTH01 £495**

Incisor hook **ANTH03 £18**

Bridge hook **ANTH02 £18**


Ceram Wire tip inc 95mm & 150mm wires **ANTH05 £60**

Simple Wire tip inc 50mm & 60mm wires **ANTH07 £46**

4 wires 2 x 50mm 2 x 60 **ANTH06 £46**

4 wires 2 x 95, 2 x 150mm **ANTH08 £46**

Adaptor 6949 **ANTH04 £34**



MISCELLANEOUS EQUIPMENT

AUTOCLAVES UK MANUFACTURE



NEW Advance Pro 16L vacuum and non vacuum cycles. Ten minutes shorter cycle than previous model. Integrated data logger. New flexi rack system can take longer, wider objects. RRP £4727

AUTO08 £3965



Alpha 16L non vacuum with 6 tray capacity. Single use water system supplied with 3 trays choice of data logger/printer. 1 yr warranty RRP 2472

ALPHA01 £2420



Classic range
21000 standard with basket. Sterilising temperature 126C time 11 mins **AUTO01 £828**
21000 gauged 134C 11 mins **AUTO19 £1112.**
Other models available

WATER DISTILLER



Water Distiller kit Includes distiller cleaner and pack of carbon sachets (not required for autoclave water) Produces 4L in 5-6 hours.

WDT01 Distiller **£169.00**

WDT02 Cleaner **£11.99**

WDT03 Storage Jar **£15.00**

WTD04 Charcoal sachets **£10.90**

HANDPIECE MAINTENANCE



Kavo Quattrocare Plus. 1 minute to clean 4 handpieces. Comes complete with 3 intramatic (E type) couplings, 1 Multiflex coupling distiller and a can of Quattrocare spray rrp £2092

Please call for our latest price and offer



iCareC2 To clean 4 handpieces. Air line required. Comes complete with 2 E type couplings, 2 M4 couplings and 1Ltr iCare oil. 500ml tank capacity. Quick connectors available for other brands **Please call for our latest price and offer**

ALGINATE MIXER



Automatic alginate mixer
Digital display touch panel control, 1-20 seconds timer can store multiple settings at the touch of a button. Patented dual start system ensures smooth start, protects gears and prolongs machine life RRP £850.00
MOTMX300 £695.00

X-RAY VIEWER

Budget X ray Viewer. View intra ora or panoramic wall mount or free standing

VIEW02 £85.00



PULP TESTER

Pulp tester ascertains the vitality alive or dead of the tooth.

DENAP01 £105.00



X RAY IMAGING SCANNER

PSPiX²



Elegant and Compact the PSPiX 2 can be installed for fast scanning directly in a surgery or as multi-user version on up to 10 workstations. High resolution editable images produced within seconds and the plate is erased ready for re use.

Includes: 4 Imaging Plates (2 x size 1, 2 x size 2) 200 bag & cover for imaging plate (10 x size 1, 100 x size 2) an imaging plate storage box, 1 stylus for touchscreen, 10 samples of Soprowipes, 1 power supply 100-240 VAC, 1 Ethernet cable & Sopro imaging software RRP £6484 our price **£4573**

PSPiX²



PSPiX 1st generation



SIZE	PSPiX2	PSPiX 1st gen
0 22 x 35mm	990215	700178
1 24 x 40mm	990216	700179
3 31 x 41mm	990217	700180
4 27 x 54mm	990218	700181
Price per pack	£60 (2)	£275 (6)

PSPiX2 bag and cover 300 pieces sizes 0-3 **£60**
PSPiX 1&2 hygiene bag 250 pieces sizes 0-3 **£49**

INSTRUMENTS

INSTRUMENT KITS



Kit1, 3 Tray Setting contains: 3 Aluminium trays, 1 pk Liners, 3 Probes, 3 Tweezers, 3 Dycal Applicators, 6 Excavators, 1 box Mirror Heads, 3 Ball Burnishers, 1 Spatula, 3 Amalgam Packers, 3 Flat Plastics, 3 Carvers. **£99.50**



Kit2, Forceps Kit contains: 1 each of 308 Upper Anterior Roots, 310 Upper Anterior 314 Upper Premolar, 326 Upper Roots Bayonet, 329 Right Upper Molars, 330 Left Upper Molars, 349 Lower Incisors, Canines & Premolars, 354 Lower Molars 361 Lower Roots. **£110.00**



Kit3 includes all items in **Kit 1 & 2** above Plus: Wide & Narrow Matrix Bands, wide & narrow Holders, Non & Self Aspirating Cartridge Syringes, Flat & Pencil Bur Brushes, Kit of 3 Elevators, Curved Scissors & Straight Scissors **£267.50**



Hygienist Kit with 4 double ended scalers H6/7, McFarlane 2/3, McFarlane 4/5, 204S, & 1 dental mirror with handle **KIT-HY £20.95**

TRAYS & RACKS




Large tray rack **DMNTR01 £59.00**
Small tray rack **DMNTR02 £49.00**

Aluminium perforated trays
DMPT__ Large tray **£5.50**
DMPTL__ Large tray lid **£5.50**
DMMP__ Small tray **£4.50**
DMMPLO1 Small tray lid **£4.75**
Large 28 x 18cm small 18 x 14cm



Aluminium non perforated trays
DMNP__ Large tray **£5.00**
DMNPLO1 Large tray lid **£5.75**
DMMP__ Small tray **£4.35**
DMMPLO1 Small tray lid **£4.50**

Waxed tray liners
Pack 250 liners
ILT001 £6.25



Tray inserts
DMALTR01 Large
DMALTR02 Small
£3.25



ARTICULATING FORCEPS



ARTF01 £7.30


AMALGAM £2.90




BDS99 (Blacks)



BDS100 (Blacks)



BDS133 (CUX5)



BDS139 (153)



AMALGAM (cont) £2.90

BDS140 (153A)



BDS141 (154)



BUR BRUSHES



Bur Brush Flat **BBF01 £1.95**



Bur Brush Pencil steel **BBS02 £3.95**

BURNISHERS £2.90



BDS119 (18)



BDS146 (155 Hopson)



BDS167 (1am-pat)



BDS168 (1Lam-pat)



CARVERS £2.90



BDS155 (Wards 1)



BDS156 (Wards 2)



BDS157 (Hollenbach)



BDS433 (Le Cron)



DYCAL APPLICATOR



DYAP01 £3.15

Equipping a new surgery? - save money with our popular instrument kits

INSTRUMENTS

FORCEPS £12.95



BDS 308 (1)
Upper Anterior Roots



BDS 310 (29)
Upper Anterior



BDS 311 (29s)
Upper Anterior



BDS 314 (7)
Upper Pre Molar Roots



BDS 319 (76)
Upper Roots



BDS 321 (76N)
Upper Roots



BDS 322 (108)
Upper Roots



BDS 324 (147)
Upper Roots



BDS 326 (52)
Upper Roots Bayonet



BDS 329 (17R)
Upper Molar Right



BDS 330 (18L)
Upper Molar Left



BDS 335 (94R)
Upper Molar Right



BDS 336 (95L)
Upper Molar Left



BDS 337 (96R)
Upper Molars Right



BDS 338 (97L)
Upper Molars Left



BDS 340 (67)
Upper Wisdom & Roots



BDS 342 (114)
Lower Centrals & Roots



BDS 345 (144) Balding
Lower Centrals & Roots



BDS 349 (75)
Lower Incisors & Canines



BDS 351 (22)
Lower Molars



BDS 354 (73)
Lower Molars



BDS 355 (73S)
Lower Molars



BDS 359 (87) Cowhorn
Lower Wisdom Molars



BDS 361 (74)
Lower Roots



BDS 364 (74N)
Lower Roots



BDS 365 (74N) Eagle Beaks
Lower Molars & Roots



BDS 366 (37)
Child's Upper Incisor Canines



BDS 367 (39)
Child's Upper Molars



BDS 368 (123)
Child's Lower Teeth & Roots



BDS 369 (138)
Child's Upper Teeth & Roots



BDS 370 (158)
Child's Upper Molars



BDS 372 (159)
Child's Upper Pre Molars



BDS 373 (160) Hawks
Child's Lower Molars



BDS 374 (110)
Child's Lower Molars



BDS 375 (162)
Child's Lower Roots Centrals



RONGEUR bone forcep

CHILDS FORCEPS £12.95



BDS 366 (37)
Child's Upper Incisor Canines



BDS 367 (39)
Child's Upper Molars



BDS 368 (123)
Child's Lower Teeth & Roots



BDS 369 (138)
Child's Upper Teeth & Roots



BDS 370 (158)
Child's Upper Molars



BDS 372 (159)
Child's Upper Pre Molars



BDS 373 (160) Hawks
Child's Lower Molars



BDS 374 (110)
Child's Lower Molars



BDS 375 (162)
Child's Lower Roots Centrals

INSTRUMENTS

Spend >£40 on BDS instruments receive 15% discount >£100 15% >£200 20%

ELEVATORS £7.90



BDS406 Cryer (L) **BDS407** Cryer (R)



Couplands 1
BDS413



Couplands 2
BDS414



Couplands 3
BDS415



Hospital Patt. **BDS393** (R) **BDS394** (L)



Lindo Levien **BDS384** (Lg) **BDS385** (Med)



Warwick James

BDS381 (R) **BDS382** (L) **BDS383** (Str)



Pereosteal
BDS411



Pereosteal
BDS412



EXCAVATORS £2.90



BDS38



BDS39



BDS40



EXCAVATORS (cont)



BDS41



BDS42



BDS48



BDS49



MATRIX BANDS



Siqueland holder £7.10
MXBH01 Wide **MXBH02** Narrow
MXB01 12 Wide Bands
MXB02 12 Narrow £0.99



Tofflemire
MXBH03 Holder £6.95
MXB03 PD assorted bands £1.75
MXB04 PD No.1 Bands £1.75

LUXATION £22.95

Stainless steel light weight <60grams



Straight sizes: **LX1S** 1mm, **LX2S** 2mm
LX3S 3mm, **LX5S** 5mm



Curved sizes 3mm, 5mm



Inverted Curve 3mm



Contra Angle 3mm
KIT of 4: 3mm & 5mm
straight and 3mm & 5mm
curved **LUX001** £85

KIT of 7 inc 2mm straight,
inverted curve & contra
angle **LUX002** £145



MIRRORS

BDS4MIRR No.
4 Plain box 12
£5.20



Mirror handle **MRH01** 0.99



MODELLING



BDS170 Osted Mitchel £2.90

ORTHODONTIC PLIERS

Long Nose Adams
BDSOP01 £12.95



Short Nose Adams
BDSOP02 £12.95



PLASTICS £2.90



BDS106
(1)



BDS113
(6)



BDS114
(6SH)



BDS117
(12 Ladmore)



BDS121
(21)



BDS123
(166)



BDS124
(176)



BDS125
(179)



BDS147
(156)



BDS148
(Smiths)



BDS154
(207/208)



INSTRUMENTS

PROBES £1.60 Unless stated

CPITN-C
£4.75



BDS04
14 Williams
£3.15



BDS06
(2)



BDS08
(6)



BDS10
(8)



BDS11
(9)



BDS21
(54)



BDS22
(Holst)



SCALERS £2.90

SINGLE ENDED

BDS 222
(Williams)



BDS 245
(Jaquette 1)



BDS 246
(Jaquette 2)



BDS 247
(Jaquette 3)



BDS 248
(4C Periodontal Hoe)



BDS 249
(7C Periodontal Hoe)



BDS 251
(9C Periodontal Hoe)



BDS 255
(H6)



BDS 256
(H7)



BDS 278
(Towner U15)



**DOUBLE ENDED WITH
LARGER STEM £4.15**

BDS 204S



BDS H6/7



BDS MF2/3



BDS MF4/5
DOUBLE ENDED £2.90

BDS 265
(A6F)



BDS 266
(Cumine 152)



BDS 268
(Shepards 197)



SCISSORS £2.90

STRAIGHT
BDS416 1 GUM
CURVED
BDS417 2 GUM



STRAIGHT
BDS418 4 GUM
CURVED
BDS419 5 GUM



STRAIGHT BDS421
CROWN BEE BEE
CURVED BDS420



STUDENT
BDS422



DRESSING
BDS423



SURGICAL STRAIGHT
BDS424S
CURVED
BDS424C



SPATULAS £2.90



ITEM 81
BDS81



20 CAPS
BDS82

SPATULAS (cont) £2.90



ALGINATE
BDS78



3 WESTERN
BDS79

SPENCER WELLS £6.00

SPENCER WELLS
STRAIGHT
SPEN02



SPENCER WELLS
CURVED
SPEN01



THYMOSIN

ITEM 2
BDS149 £2.90



ITEM 3
BDS150 £2.90



TWEEZERS

£2.90



PERRY
SERRATED
BDS425



COLLEGE
BDS426



PLUGGING
BDS427



GUTTMAN
SERRATED
BDS428



SELF LOCK
SELF01



**RAT TOOTH
TISSUE FORCEPS
RAT01 £2.90**

WAX KNIVES £2.90



BDS430 SMALL KNIFE



BDS431 LARGE KNIFE

**WILLIS BITE
GAUGE**



**WIRE
CUTTERS**



CONSUMABLES

ALGINATES

Kromopan 450g high precision and elastic impression chromatic phase indicator
KA01 bag **£6.99**
Box rate £6.39 each



Empress 454g very accurate detailed reproduction, elastic fast set 130 secs.
 total time **EA02 £5.45**
10 bag rate £5.07



Hydrogum5
 5 day high stability, lilac colour with fruit berry flavour rapid 45 sec.
 set time **HYDR05 £8.65**



Pegasus trays

Boxes of 25
 1 Large upper
 2 Large lower
 3 Medium upper
 4 Medium lower
 5 Small upper
 6 Small lower
 7 Partial up Left
 8 Partial up right
 9 Anterior
IMPRO_ £5.45



Unodent trays

Boxes of 25 **SITXX £5.50**
 3 Upper deep eden.
 4 Lower deep eden.
 5 Upper shallow eden.
 6 Lower shallow eden.
 11 Upper dentate lg. 12 Lower dentate lg.
 13 Upper dentate av. 14 Lower dentate av.
 15 Upper dentate ort. 16 Lower dentate ort.
 17 Upper dentate sm. 18 Lower dentate sm.



Dentsply polytrays Boxes of 25 **£19.24**
 Lower:

POLYTL0 sm. green **POLYTL1** med. blue
POLYTL2 lg. yellow

Upper:

POLYTU0 sm. green **POLYTU1** med. blue
POLYTU2 lg. yellow

Alginate mixing bowls & spatulas

ALB01 bowl **£2.95**
ALS01 spatula **£2.95**



IMPRESSION MATERIALS



Bonasil A+ wash Low viscosity
 A silicon impression material, thixotropic 2 x 50ml + tips. Very high tear strength
BONAW01 regular **£17.95**
BONAW02 fast **£17.95**



Bonasil A+ putty High viscosity putty
 400g base 400g catalyst, regular (work 2 min set 2.5 min) and fast (work 1 min set 2 min)
BONAP01 regular set **£29.95**
BONAP02 fast set **£29.95**

Tray tack adhesive
 Spray on 200ml can
TTS01 £8.65
 Paint on 30ml bottle
TTP01 £6.45



Pegasus orange spray 200ml solvent
 spray aids cleaning and removal of alginate and adhesive from re-usable trays
PEGO01 £5.75



ETCHANTS



BDSI ETCHANT GEL UK manufacture
 40% aqueous orthophosphoric acid green in 50ml syringe for dispensing to 10 applicator syringes + 60 tips
ETCH04 50ml **£22.50 BEST VALUE**



PEGASUS ETCHANT GEL 35% phosphoric acid green blue gel in 12gm syringe 20 tips **ETCH02A £8.50**
 3 rate **£8.00** 5 rate **£7.50**

AMALGAM



ULTRACAPS+ Economical high strength admix amalgam Ag 50% Sn 30% Cu 20% in regular and fast set

Boxes of 50 capsules **ULTC1** 1sp **£27.56**
ULTC2 2sp **£35.64**, **ULTC3** 3sp **£41.41**
ULTC5 5sp **£67.79**.
 Tubs of 500 capsules **BULTC1** 1sp **£232.35**
BULTC2 2sp **£304.16** **BULTC3** 3sp **£380.12**
BULTC5 5sp **£546.81**

GS80 Economical high strength admix amalgam Ag 40% Sn 31% Cu 29% in regular and fast set

Boxes of 50 capsules **SPIL01** 1sp **£31.85**
SPIL02 2sp **£39.88**, **SPIL03** 3sp **£47.56**
SPIL05 5sp **£71.22**
 Tubs of 500 capsules **BSPLI1** 1sp **£279.00**
BSPLI2 2sp **£342.09** **BSPLI3** 3sp **£396.13**
BSPLI5 5sp **£627.70**

GS80 Spherical admix amalgam
 Ag 56% Sn 29% Cu 15%

Boxes of 50 capsules **SSPIL1** 1sp **£33.47**
SSPIL2 2sp **£41.39**, **SSPIL3** 3sp **£49.49**
SSPIL5 5sp **£74.26**
 Tubs of 500 capsules **SBSPLI1** 1sp **£303.98**
SBSPLI2 2sp **£359.29** **SBSPLI3** 3sp **£444.88**
SBSPLI5 5sp **£658.87**

LOGIC+ High silver high strength spherical amalgam Ag 60% Sn 28% Cu 12%

Boxes of 50 capsules **LSPIL1** 1sp **£46.75**
LSPIL2 2sp **£58.65** **LSPIL3** 3sp **£68.40**
LSPIL5 5sp **£101.45**
 Tubs of 500 capsules **BLSPIL1** 1sp **£372.61**
BLSPIL2 2sp **£473.03** **BLSPIL3** 3sp **£559.56**
BLSPIL5 5sp **£848.08**

PERMITE High silver high strength admix amalgam Ag 56% Sn 28% Cu 15%

Boxes of 50 capsules **PSPIL1** 1sp **£46.75**
PSPIL2 2sp **£58.65** **PSPIL3** 3sp **£68.40**
PSPIL5 5sp **£101.45**
 Tubs of 500 capsules **BPSPIL1** 1sp **£372.61**
BPSPIL2 2sp **£473.03** **BPSPIL3** 3sp **£559.56**
BPSPIL5 5sp **£848.08**

BUY ANY 3 SDI 500 TUBS AND GET AN SDI ULTRAMAT S AMALGAMATOR FREE

ETCHANT TIPS



Metal Etch Tips MET01 pk 50 **£8.95**

CONSUMABLES

COMPOSITES



BDI Composite in compules or syringes
UK manufacture, fluoride release, excellent radiopacity for I,II,III,IV,V

Compules 20 x 0.25g **£23.35 5 rate £22.89**
4 shades **COM248X2A1** Light,
COM248X2A2 Light yellow, **COM248X2A3**
Universal, **COM248X2A3.5** Dark Yellow

BEST SELLER - LOW PRICE




4g syringes **£10.22 5 rate £9.52**
4 shades **COM248X3A1** Light,
COM248X3A2 Light yellow, **COM248X3A3**
Universal, **COM248X3A3.5** Dark Yellow



BRIGHT LIGHT FLOWABLE COMPOSITE
A light cured microhybrid low viscosity flowable composite in 2gm syringe, 4 shades **FCA1, A2, A3, A3.5** **£4.95**



BRIGHT LIGHT PASTE/PASTE COMPOSITE
Self cured kit 14g base paste, 14g catalyst, pads & spatulas **PASTE02** **£19.95**



Universal nanohybrid composite. Simplifies restorations, containing four simple shades that equate to the Vita shade guide.

Compleat Kit 20x4 0.25g **£172.37**
Refills 20 x0.25g **AUEA1-3.5** **£45.00**
Syringe Kit 4x1 4g **£145.05**
Refills 1 x4g **£37.95**

GLASS IONOMERS



VERY POPULAR - GREAT VALUE

CERAMFIL B GI Kit for class III & V
UK manufactured kit 5 x 10gm pots in 7 shades. **CERK1** Ivory, 2 Light Yellow, 3 Universal, 4 Light Grey, 5 Yellow, 6 Light Brown, 7 Dark Yellow. Combination kit contains 1,2,4,5 & 7
Per Kit **£52.50 2 rate £51.00**
5 rate £45.50 Single Pots 12.95



CERAMFIL MOLAR A high strength condensable microglass ionomer perfect for small **posterior** fillings 25g powder 15ml liquid **CRMM01** **£26.95**
BEST SELLER



Ceramsave B Aquaset Fast setting GI for restoration of **deciduous** teeth, fluoride release, radiopaque. **CMS01** 30gm powder & accs. **£21.50**
2 rate £19.95



Ceram II VLC Light cured GI class III & V for **anterior** teeth, fluoride release, high translucency. **CERA01** 10gm powder & accs. **£24.95**



Ceramfil A.R.T. GI for class III & V especially developed for the **Altraumatic Restorative Technique**. 30gm powder, 9ml varnish & accessories
CART01 **£25.85**



Ceramcore B Aquaset Fast setting GI for fabrication of cores and tooth build up, fluoride release, radiopaque. Easy to mix & trim as part of crown prep during same visit
30gm powder & accs. **CNC01** **£19.95**

Ceramcore B Silver Fast setting GI high compressive strength, excellent abrasion due to silver content, bland to pulp contours & carves easily. **SRI01** 30gm powder & accessories **£29.50**

RIVA



Riva Self Cure High strength GI in regular & fast set, fluoride release. Boxes of 50 capsules in 7 shades Regular set: **RIVARSA1, A2, A3, A3.5, A4, B2, B3** **£73.30**
Fast set **RIVAFS**



Riva Light Cure Resin reinforced GI, high fluoride release. Boxes of 50 capsules in 11 shades **RIVALCA1, A2, A3, A3.5, A4, B2, B3, B4, C2, C4, b** **£81.50**



Riva Silver Silver reinforced GI mainly for Class II & III, deciduous and geriatric teeth and core build ups. Boxes of 50 capsules **RIVASIL01** **£102.15**



BDI composite made in the UK is our best selling composite

CONSUMABLES

CEMENTS & LINERS



Orthocem B Unique GI bonding adhesive for use in cementation of bands. Fluoride release, no acid etch, translucent, bonds metal to tooth, bio compatible
OCM01 30gm powder & accessories £22.50



Ceramchem Aquaset Fast set permanent luting cement for cementation of crowns fluorides release, high compressive strength
CMC01 30gm powder & accessories 26.95



Carbchem Aquaset Fast setting polycarboxylate luting cement. Can be used for durable temporary filling
CBC01 60gm powder & accessories £19.95



Ceram I New GI luting cement powder and liquid formulation. Ideal for cementation and also lining dentine wall of cavities for composite restoration
CERIO1 25g powder 15ml liquid £23.95

BDSI

... for friendly, helpful service



seT all in 1 luting cement 5 shades A1, A2, OA3, white opaque & translucent
SET__ (shade) box of 50 capsules £139
SETpp 2 x 7g refill syringes £76.15



Riva Luting A high fluoride releasing GI cement for cementation
RIVALC box 50 capsules £57.75



Riva Protect A high fluoride releasing GI cement for surface protection, sealing and caries stabilisation
RIVAPRO01 (bleach shade) box 50 £81.75
RIVAPRO02 (pink shade) box 50 £81.75



Riva Luting Plus A high fluoride releasing GI luting cement for final cementation of metal PFM and resin crowns, radiopaque, universal light yellow shade
RIVALPC box 50 capsules £71.44



BRIGHT Cavity Liner Glass reinforced light cure single component cavity liner, fluoride release, radiopaque to enforce the floor and walls of the cavity prior to filling
CL01 2gm syringe £9.95



BDSI Zinc Oxide UK manufacture reinforced Eugenol filling material
KLZ03 40g powder £7.99
KLZ04 15ml liquid £7.99



KALZINOL Resin bonded zinc oxide eugenol cement
KLZ01 40g powder £21.45
KLZ02 15ml liquid £18.95



DYCAL Radiopaque calcium hydroxide cavity lining material
DYC01 13g base 11g catalyst £23.79

ADHESIVES



BDSI ONE Step Bond UK Manufactured seventh generation light cure self etching one step bonding adhesive.
BOND01 6ml + accessories £29.95



go! A single component light cured self etching bond
GO01 1 x 5ml bottle £66.70
GO02 50 x 0.1ml single unidose kit £69.95

CONSUMABLES

POLISHING



Prophy cups latch 100

PCL01 £10.50

5 bag rate £9.50 10 bag rate £8.50



Prophy cups snap on 100

PCSN01 £7.50

5 bag rate £7.00 10 bag rate £6.50



Prophy cups screw in 100

PCS01 £9.50

5 bag rate £9.00 10 bag rate £8.50



Bristle latch 100 **BBL01 £17**

Nylon latch 100 **BBL01N £17**



Bristle screw in 100 **BBS01 £17**

Nylon screw in 100 **BBS01N £17**

Bristle snap on 100 **BBSN01 £18.00**

Prophy Paste

mint flavour

medium grain,

excellent

quality 250gm jar

PROPH01 £4.95

BUY 5 at £4.25 each



Acrylic Trimmers ISO size 060



Bud HP Cross Cut

Ref 85

ACR02 £11.85



Pear HP Cross Cut

Ref 87

ACR06 £11.85



Bud HP Plain Cut

Ref 75

ACR04 £11.35



Pear HP Plain Cut

Ref 77

ACR07 £11.35

BURS

Group 1 bur discounts:

Buy 20 get 2 FREE = 86p each

Buy 30 get 4 FREE = 84p each

BUY 50 get 7 FREE = 83p each

Diamond Group 1 95p each



501 502 503 504 512 520 521 522 525 526 527 530 531 535 536 537 538 540 541 542



543 544 545 546 550 551 552 554 555 556 557 558 559 560 561 562 563 564 582

Diamond Group 2 £1.35




565 566 567 607 608 613 615 622 624 626 636 653 654 700 708 709 721 724 730 731 745

Diamond Group 3 £1.69



629 631 632 637 640

Steel Burs latch fitting in packs of 6



Flat Fissure sizes: **FFSB00**

1/2, 1, 2, 3, 4, 5, 6, 7 **£2.40**

Round sizes: **RSB00**

1/2, 1, 2, 3, 4, 5, 6, 7, 8 **£2.40**

9, 10, 11, 12 **£4.50**

Inverted Cone sizes: **ICSB00**

1/2, 1, 2, 3, 4, 5, 6 **£2.40**

Diamond Group 4 £2.25



513 514 515 592 0.3

Tungsten Carbide Burs Pack of 5 £5.75

ISO Ø 1/10mm	010	014	016	ISO Ø 1/10mm	010	012
Head L mm	1.2	1.4	1.6	Head L mm	4.4	4.5
US No.	35	37	38	US No.	1157	1158
TC_ _ (FG)	35	37	38	TC_ _ (FG)	1157	1158

ISO Ø 1/10mm	009	012
Head L mm	4.0	4.5
US No.	56	58
TC_ _ (FG)	56	58

ISO Ø 1/10mm	010	012
Head L mm	4.0	4.5
US No.	1557	1558
TC_ _ (FG)	1557	1558

ISO Ø 1/10mm	012	008	010
Head L mm	2.0	1.0	1.5
US No.	332	330	331
TC_ _ (FG)	332	330	331

ISO Ø 1/10mm	012	016	021
Head L mm	4.5	5.0	5.5
US No.	701	702	703
TC_ _ (FG)	701	702	703
TC_ _ (RA)	701	702	-
TC_ _ (HP)	701	702	-

ISO Ø 1/10mm	009	010	012
Head L mm	4.0	4.5	4.5
US No.	556	557	558
TC_ _ (FG)	556	557	558
TC_ _ (RA)	-	557	558
TC_ _ (HP)	-	557	-

ISO Ø 1/10mm	008	010	012	014	016	018
Head L mm	X	X	X	X	X	X
US No.	1	2	3	4	5	6
TC_ _ (FG)	01	02	03	04	05	06
TC_ _ (RA)	01	02	03	04	05	06

CONSUMABLES

ENDO FILES



NiTi Greater Taper Files 4% or 6% taper
25mm long, packs of 6 sizes **15** white, **20** yell,
25 red, **30** blue, **35** green, **40** black & Asst
15-40 speed 150-350 rpm
NiTi4%__ **NiTi6%**__ **£12.99**



NiTi Multi Taper U Files Shaping &
Finishing files 25mm long, (19mm for MSX)
Packs of 6 shaping **MSX**, **MS1** purple, **MS2**
white finishing **MF1** yell, **MF2** red, **MF3** blue,
Asst 1 of each MSX to MF3
NiTi__ **£17.99**

Great value very popular



Wave+ Gold One file system, Asymmetric
parallelogram cross section, 1-2 cutting
points acting on root canal wall. Progressive
taper. Crown down technique, reciprocating.
Packs of 6, 25mm long **20** small, **25** primary,
35 medium, **40** large
NiTi__ **£19.99**



H Files 25mm long sizes 10-55, 60, 70, 80
HFI__ **£3.75** **10 rate £3.25**



K Files 25 & 21mm sizes 10-55, 60, 70, 80
KFI__ **£3.75** **10 rate £3.25**



Peeso Reamers 1,2,3,4,5,6 **PR**__ **£7.25**



Gates Glidden 1,2,3,4,5,6 **GG**__ **£7.95**

Hand Reamers 25mm **HRE**__ **£3.75** **Finger Pluggers** 25mm
Asst **FP15-40** **£3.75**

Paste Carriers Pk4 25, 30, 35, 40
Asst 25-40 **PF**__ **£3.75** **Finger Spreaders**
25mm, 15, 20, 25, 30
Asst 15-40 **FS**__ **£5.75**

GP POINTS

GP Points colour
coded hand rolled
120 per pack sizes
10-90p **GP**__
£5.25



GP Taper Points
.06 and .04
taper 60 per pack
sizes 20-40 **GP6%**__
£6.50 **GP4%**__ **£6.50**



GP W-one points
Small, primary, med. &
large 60 per pack
GPW1__ **£9.95**



PAPER POINTS

Paper Points colour
coded 100 per pack
sizes 10-90pm **PAP**__
£3.99 **3+ rate £3.51**



Paper Taper Points
.06 and .04 taper
100 per pack
sizes 20-40 **PAP6%**__
£4.99 **PAP4%**__ **£4.99**



Paper Points W-one
Small, primary, med. &
large 100 per pack
PAPW1__ **£8.95**



MIXING ITEMS



Mixing pads
Size A 2.5" x 1.5" Pack 5 **MPA01** **£3.85**
Size B 3" x 2.5" Pack 5 **MPB01** **£7.95**
Size C 6" x 3.5" Pack 5 **MPC01** **£8.95**



Glass slabs
Plain 4" x 4" x 1/2" **GSPP01** **£3.90**
Plain 6" x 3" x 5/8" **GSPP02** **£10.25**
Frosted 4" x 4" x 1/2" **GSF01** **£3.90**

Plastic mixing wells
12 x 12 wells
MXW01 **£3.75**



MIXING ITEMS



Dappens pots
Glass: Clear **GDP02**, Blue **GDP03**,
Amber **GD4** Green **GDP05** **£2.99**
Plastic: packs of 200, 50 of each **£7.95**

SDI points micro
applicators, boxes 400
POIN01 medium.
POIN02 assorted
POIN03 fine
POIN04 superfine
£19.75



Economical micro
applicators, boxes
400 regular
MAPU01 **£17.50**



Mixing spatulas pack
100 **MXS01** **£4.55**



MOUTHWASH & CUPS



BDSI mouthwash concentrate 100ml
Spearmint flavour **MWC01** **£9.99**
6 rate £8.99



BDSI mouthwash tablets 1000 per
jar in Thymol, Mint Spearmint and
Mixed **MWT01** **£5.95**



Disposable Cups 7oz squat 2000
per box **DC01** **£23.95**
OFFER Buy with 2000 mouthwash
tablets DC02 **£32.99**

CONSUMABLES

DISINFECTION



Desderman pure gel advanced antimicrobial hand disinfection kills 99.99% of viruses and bacteria leaving skin soft 500ml **DESDER01** hand rub **£5.79**
DESDER03 pure gel **£5.79**



Mikroqid wipes 200 wipes. Hard surface wipes effective against bacteria, fungi and viruses inc HIV and MRSA
MIKR01 tub alcohol free **£7.65**
MIKR02 refill alcohol free **£7.30**
MIKR06 tub alcohol **£5.99**
MIKR07 refill alcohol **£6.65**



Gigasonic
High performance general purpose cleaning and disinfection solution for ultrasonic baths 2Ltr. Dilute 3% in use
GIGA02 (rrp 60.63) **£47.95**



Classic Hard surface disinfectant spray economy refill **SURFD01** **£17.90**



AspiSept
Aldehyde free concentrated solution for cleaning and disinfecting dental aspirating units. Dilute 2% with tap water
ASPI01 **£16.95**



Gigazyme plus Manual high level disinfection of dental instruments prior to sterilisation. Can also be used in ultrasonic baths. 1Ltr dilutes to 200Lts
GIGA03 **£22.16**

SURGERY ITEMS

Examination Gloves



Our gloves are certified to BS EN455 parts 1,2 & 3. Part 3 ensures they are fully compliant

Blue Nitrile powder free non latex
200 box Xsmall **MINI01**, small **MINI02**, medium **MINI03**, large **MINI04**, Xlarge **MINI05** **£7.99** **BUY 10 BOXES at £5.99**
BUY 20 BOXES at £4.99 each

Face Masks



Tie on mask 50 FM02 Loop 50 FM01
Contact us for price and availability



Blue velcro bibs Patient fastening velcro bibs **BVB01** **£2.05** **10 rate £1.85**



Disposable bibs
Two ply polythene backed, absorbent tissue, rectangular bibs for use with bib chains **DISAPR01** box 500 **£17.50**



Plastic aprons 27" x 42" packs of 100
DA01 white **£7.95** **DA02** blue **£7.95**
5 rate £6.95



Mini centre Feed Rolls 2 ply white
12 pack 195 x 400mm **MCF01**
£22.00 **2 rate £18.00**



C Fold 2 ply white luxury towels
210x230mm **HT-01** 2400 towels
£18.00 **BUY 2 at £16.00 each.**



Cotton Wool Rolls 1000 per box
CWR03 size 1 **CWR01** size 2
£5.95 **5 rate £4.95**



Gauze Napkins 12.5 x 12.5
500 per box **GAUZ01** **£6.95**
3 rate £6.25



Throat Pack 10 mtr roll
THR01 **£4.40**



Paper aprons 50 x 80cm with perforated cut out for the neck. Tie fastening **BIBR01** **£8.95**

SURGERY ITEMS



Facial Tissues High quality 200 x 210
100 per box, 24 boxes per case
PROF01 £26.97 **2 rate** £21.95



Autoclave Pouches boxes 200
Granton Medical UK manufacture
ACP05 90 x 230 £4.95 **new low price**
ACP02 90 x 200 £5.95
ACP10 135 x 255 £9.75



Cyberpouch self seal boxes 200
ACP08 83 x 159 £7.50



Saliva ejectors transparent flexible
clear plastic. Metal strip enables
ejector to bend. Pack of 100.
Smooth edges and comfortable
13 cm long. **SEC01** £1.95



Disposable aspirator tips Pack of 100.
Vented **MAT01** £2.30

Barrier covers



Pen type for curing lights and air/water
syringes 64mm x 255mm pack 500
AWSC01 £9.95



Tubing covers pack 250
Large 457 x 53mm **TUBC01** £9.95
Small 457 x 44mm **TUBC02** £9.95



Head rest cover 285mm x 255mm
pack 250 **HRDBC01** £6.50



Half chair cover 840 x 916mm
pack 125 **HCC01** £19.95



X-ray sensor cover 222 x 41mm
pack 500 **XRSC01** £5.95



Digital X-ray barrier covers
Size 0 20x30 (100) **UNOXSC000** £11.50
Size 1 20x40 (100) **UNOXSC001** £14.50
Size 2 30x40 (300) **UNOXSC002** £30.00
Size 3 57x54 (100) **UNOXSC003** £16.50
Size 4 54x75 (100) **UNOXSC004** £24.00

For PSPIX plates and bags
please see page 14

SUNDRY ITEMS



Wraparound safety glasses. With side
shields 2 position adjustable S/M/L
SGA01 £4.95



Hogies + Eyeguards Stylish with
adjustable nosepiece. Ladies black blue,
pink. Gents black, blue. Childs blue
BLHOH01 £29.50



Amalgam carriers **AMC01** straight
AMC02 45 degree **AMC03** 90 degree
£5.90 each **6 rate** £5.70

3 IN 1 SYRINGES
DCI type air & water
syringe with 2 syringe
tips **3IN1S01** £39.95
Pack of 5 metal tips
3IN1ST01 £9.95



Disposable 3 in 1
tips **DISTIP** box 200
£16.50 **5 rate** 14.95



Cartridge syringes 2.2ml
CSA01 non aspirating
£9.95 **5 rate** £8.95
CSA02 self aspirating
£12.95 **5 rate** 11.95



Tubing Adaptors
ADAP01 connects Borden 2/3 hole
handpiece to midwest tubing
ADAPT02 connects Mid West 4/5 hole
handpiece to borden tubing
£39.95 each

TOOTH WHITENING

BOUTIQUE WHITENING

Boutique Whitening is a premium quality whitening gel that delivers superior results with a noticeable improvement in just a few days.

Boutique By Day, whitens with only 1 hour wear time per day. Alternatively Boutique By Night whiten your patients teeth as they sleep



Boutique 4 x 3gm syringe kits **£22.50**
BOUT6KIT Hydrogen Peroxide 6%
BOUT10KIT Carbamide Peroxide 10%
BOUT16KIT Carbamide Peroxide 16%



Top up 2 syringe pillows **£10.00**
BOUT62P 2 syringes 6% HP
BOUT102P 2 syringes 10% CP
BOUT162P 2 syringes 16% CP



Voco Profluorid Varnish
 Desensitizing for hypersensitive teeth
 10ml tubes
VOPV02 Cherry **£16.11 new lower price**



Soothe desensitizing Gel
 6% Potassium Nitrate 0.1% fluoride
 gel. 4 x 1.2ml syringes
SOOTH1 **£36.22**



The high viscosity, neutral pH tooth whitening gels ensure the greatest patient comfort in a take home kit. The high water content assists in reducing sensitivity.



POLADAY



Poladay 6% Hydrogen Peroxide in 1.3g syringes.

Treatment guidelines for day
 3% 2 or 3 x 30min/day or
 1 x 60 min/day
 6% 1 or 2 x 30 min/day or
 1 x 45 min/day

Mini 4 syringe kits:
POLADM6 6% **£25.95 3@£23.99**

Standard 10 Syringe kits:
POLAD6 6% **£45.10 3@£42.99**

Bulk 50 syringe kits
POLADB6 **£165.90 3@£156.50**



BUY 3+ Kits at our discounted rates and in addition get the latest SDI special offers PLEASE CALL 01226 380400

POLANIGHT



Polanight 10% or 16% Carbamide Peroxide in 1.3g syringes

Treatment guidelines for night 10% 2 hrs to overnight 16% 90 min to overnight

Mini 4 syringe kits: **£25.95**
3@£23.99
POLANM10 10% **POLANM16** 16%

Standard 10 Syringe kits: **£45.10**
3@£42.99
POLAN10 10% kit **POLAN16** 16%

Bulk 50 syringe kits **£165.90**
3@£156.50
POLANB 10% kit **POLANB16** 16% kit



Pola Zip Lock Pouch
 a presentable way to dispense loose bulk syringes
POLAP02 **£18.75**



POLA OFFICE



Polanight 10% or 16% Carbamide Peroxide. No bleaching light required, the built in desensitiser & high pH ensures maximum comfort. Simple delivery with dual barrel syringes. In 1 and 3 patient kits

1 Patient kit: 1 x 2.8ml syringe and gingival barrier **POLAOFF61** **£32.00**
 3 Patient kit **POLAOFF63** **£87.50**

STOOLS

Highly polished 5 star base with twin chrome hooded castors.
Individually assembled and stitched in the UK.



Assistant stool 9" gas lift with
foot ring **GRA11** £215
Surgeon stool no ring
GRS12 £175



Surgeon stool 9" gas lift with
foot ring **GRA14** £265
Surgeon stool no ring
GRS13 £225



Surgeon stool with centre
curved backrest **SRS08** £235



Assistants stool with 12" gas
lift backrest either left/right
position **SRA09** £270



Assistants stool 12" gas lift
extra lever for back rest free
float action, lockable position
SRA07 £325
Surgeon stool no foot ring 9"
gas lift **SRS06** £295



Surgeon saddle stool,
ergonomically designed to
retain correct posture, tilt
mechanism 9" gas lift
GMS10 £315



Assistant stool back support
or arm rest J arm rotates 360°
and height adjustable arm rest
fixed left/right on a ratchet
12" gas lift **SRA05** £335

Dental stools made to order. Choose your style and colour online at

www.ukdentalsupplies.com/product-category/equipment/stools

Prices are for Naugahyde colours. Belmont and ADEC also available



Unit 14 Claycliffe Business Park, Barugh Green, Barnsley, S75 1JU
Tel: 01226 380400 email: sales@bdsi-ltd.co.uk
website: ukdentalsupplies.com
Hours 8.30am -5.00pm Monday - Friday

All prices subject to the current rate of vat.
We reserve the right to correct clerical E&O.
Full terms and conditions available on request or
alternatively view online at ukdentalsupplies.com

