



Paul Tipton



Tracy Doole



NSK

TAKE THE
LEADING EDGE.

2020



Gemma O'Callaghan



Dr Nilesh R. Parmar

SPECIAL
OFFERS
INSIDE.

NSK UK Ltd
www.nsk-uk.com
0800 6341909



Dr Rhodri Thomas



INSPIRED
BY US.

NSK PRODUCTS ARE
FOR DENTISTS,
HYGIENISTS &
PRACTITIONERS
OF EVERY KIND.

Introducing Osseo 100

30% off
Only
£1,450*

Osseointegration Monitoring Device

Osseo100 RRP £2,075
Code Y1004175

Measure implant stability and osseointegration to enhance decisions about when to load the implant. Especially important when using protocols with shorter treatment time and treating higher risk patients.

Speak to your local NSK representative today to find out more



3-step procedure



R RENTAL
FROM **NSK**

Considered renting? Contact us on
0800 634 1909 or visit nskrental.co.uk

Your Local Product Specialists


**Public Sector &
Decon Manager**
Mark Beckwith
07900 246529
mbeckwith@nsk-uk.com

N.Ireland & ROI
Product Specialist / Technical Services Engineer
Jonathan Singh
UK: 07464 675158
ROI: 01695 0053
jsingh@nsk-uk.com

North West/North Wales
Andy Hargreaves
07872 824175
ahargreaves@nsk-uk.com

Contact Us

Office 4
Gateway 1000
Arlington Business Park
Whittle Way
Stevenage
Hertfordshire
SG1 2FP



Scotland
Territory Manager
Angela Glasgow
07525 911006
aglasgow@nsk-uk.com

**Midlands &
The North East**
Graham Brown
07541 864641
gbrown@nsk-uk.com

Yorkshire
Jo Dronfield
07545 786619
jdronfield@nsk-uk.com

South East
Francois Faro
07900 246822
ffaro@nsk-uk.com

London W1, WC1/2
Fiona Breitenbach
07590 443835
fbreitenbach@nsk-uk.com

South (Central)
Steve Barnes
07889 181099
sbarnes@nsk-uk.com

South West/ South Wales
Ralf Wendel
07391 015229
rwendel@nsk-uk.com



Call us on 0800 6341909 or visit www.nsk-uk.com



[nsk.uk](https://www.instagram.com/nsk_uk)



NSK UK Ltd



@NSK_UK

Are you prepped for digital? Go Electric with NSK

Accurate and clean margins are the most important factor when achieving the best digital scan. Electric driven speed increasing handpieces are the best option when looking for the control and torque you need to produce accurate margins.

What is a speed increasing handpiece?

A contra-angle that does everything a turbine (high speed) does and uses the same FG burs. BUT it has incredible 'torque' originating from the micromotor which puts it in a league of its own and puts you in total control - no more loss of power as you come into contact with the tooth. Choose your speed. Save time cutting preps.

What motor do I need to use speed increasing handpieces?

It must be an electric micromotor (NOT an air motor) as the gears inside (usually 1:5) need the electrical power to run at a high speed. These ultimate brushless electric micromotors can be retro fitted by an engineer or come as an option when buying new chairs. NSK's NLX micromotors offer a powerful 3.4Ncm of torque.

Can I run my old contra-angles (slow speed handpieces) on it?

Yes, but with a micromotor and a control panel you have the option to slow the speed right down. You can now also use handpieces with internal spray and fibre optics.

Will there be any wires showing?

A built-in version will not show any wires.

Always check compatibility to upgrade.



Options for Speed Increasing Handpieces (1:5)

	S-Max		Ti-Max
Model	M95 (non-optic)	M95L (optic)	Z95L (optic)
Lightweight titanium outer casing			X
Narrower neck/body			X
4 hole spray	X	X	X
Warranty	2 years†	2 years†	3 years
RRP	£775	£ 981	£1446



NSK NLX micromotors are available on A-dec and Belmont equipment



Built-in Kit Options

At the touch of a button, control your micromotor and/or ultrasonic scaler. All built-in kit options come with built-in kit and will need a qualified engineer to fit (additional cost).

FREE Micromotor & built-in kit (option 1) when you buy 3 x Ti-Max Z95L*

*Can include Nano95LS (as Z95L but 10% shorter) or Ti-Max Z45L

Alternative offer: Buy 4 Ti-Max Z Series contra-angles and get option 2 FREE*

Alternative offer: Buy any of the Built-in or Plug & Play options below and get a FREE nano95LS speed increasing handpiece*

- 1. RRP £2,020
Multipad (Control pad) + NLX nano (optic micromotor)
- 2. RRP £2,616
Multipad + NLX Plus (optic micromotor also with Endo features – auto reverse on/off and auto stop)
- 3. RRP £3,363
As option 1 + Varios 170 LED (built-in Ultrasonic handpiece with tips for scaling, perio, endo, restorative work)
- 4. RRP £3,957
As option 2 + Varios 170 LED
- 5. RRP £1,416
Varios 170 LUX + FREE Multipad
- 6. RRP £2,383
NLX Nano portable (Plug & Play)

NLX nano portable



2 year warranty

Get prepped for digital at www.NSKgoelectric.co.uk

Calling all BienAir MC3® (brush motor) users!



Buy an M40XS LED optic brushed motor and Ti-Max X95L for £1,200+VAT*
Buy an M40N XS non-optic brushed motor and S-Max M95 for £975+VAT*

*Terms and conditions apply. Offers only available through participating dealers. †With effect of 1st January 2019 for products sold by NSK United Kingdom Ltd and used in the UK and RoI. MC3® is a registered trademark of BienAir®.

Contra-angles

A selection of NSK Speed Increasing, 1:1 Direct Drive, Speed Reducing and Straight handpieces. Choose the features you require. The basic range of FX or EX (12 month warranty), then moving up to the Advanced range of S-Max (18 month warranty) all with internal spray and options with or without fibre optics. Move into the Premium Ti-Max series which have the added advantages of being lightweight titanium. The necks and heads become smaller as you move into X Series (2 year warranty.) and furthermore as you move into the Premium Z Series (3 year warranty).

Ti-Max Z

Buy 2, get an extra 1 FREE*



††Only external body component is titanium

- Ti-Max**

 - Lightweight titanium external body††
 - Small head 8.9mm, H13.4mm (Z series)
- Scratch resistant with DURAGRIP (Z series)
 - Option of 45° angled head (Z45L)
 - Nano range (10% shorter and lighter).

Ti-Max Z Series – Premium Contra-Angles					
Optic	Speed / Type	Model	Order code	RRP	Warranty
	1:5 increasing	Z95L	C1034001	£1,446	3 years
	1:5 increasing	nano 95LS	C1099001	£1,446	
	1:4:2 increasing and 45° head	Z45L	C1064001	£1,446	
	1:1	Z25L	C1038001	£1,203	

S-Max M

FREE UPGRADE*
from S-Max M to Ti-Max Z when you trade in a turbine **

MORE TORQUE

Buy 3, get an extra 1 FREE*



- S-Max**
- Stainless steel body
 - Superb grip
 - Powerful cutting

S-Max - Advanced Contra-Angles					
Optic	Speed / Type	Model	Order code	RRP	Warranty
	1:5 increasing	M95L	C1023001	£981	2 years†
	1:1	M25L	C1024001	£794	
	4:1 reducing	M15L	C1025001	£929	
Non-optic	1:5 increasing	M95	C1026001	£775	2 years†
	1:1	M25	C1027001	£522	
	4:1 reducing	M15	C1028001	£620	
	1:1 straight h/p	M65	H1008001	£486	

FX Series

NAC-EC THE ORIGINAL

Buy 1 FX get 1 FREE*



NAC-EC Triple pack offer £180*

FX & NAC-EC

The original NAC-EC is still commonly used. However the FX25m offers a more robust alternative, with a push button chuck and a longer neck for easier access for posteriors.

FX and EX Series – Budget Contra-Angles					
Non-optic	Speed / Type	Model	Order code	RRP	Warranty
	1:1	FX25m	C1052001	£239	12 months
	1:1	FX65m	H1014001	£239	
	1:1	NAC-EC	Y110148	£141	

*Terms and conditions apply. Offers only available through participating dealers. †With effect of 1st January 2019 for products sold by NSK United Kingdom Ltd and used in the UK and RoI. **One per instrument



Ti·Max Z



Buy 2
& get
extra 1
FREE*

Optic	Ti-Max Premium Turbines					
	Fitting	Model	Size	Order Code	RRP	Warranty
	NSK	Z900L	Std head: Ø 12.5 x 13.1mm	P1111001	£1106	3 years
	KaVo®	Z900KL		P1113001		
	Sirona®	Z900SL		P1115001		
	W&H®	Z900WL		P1121001		
	Bien-Air®	Z900BL		P1117001		
	NSK	Z800L	Min head: Ø 10.8 x 12.1mm	P1110001	£1106	
	KaVo®	Z800KL		P1112001		
	Sirona®	Z800SL		P1114001		
	W&H®	Z800WL		P1120001		
	Bien-Air®	Z800BL		P1116001		

- Our most powerful turbine range
- 94° angled head for easier access
- Lightweight titanium body**

††Only external body component is titanium

"I use a range of NSK handpieces, and in particular the NSK Z900L air turbine. The LED optics are exceptional and the Z900L takes my restorative work to a new level as it delivers a powerful and responsive performance and best of all, it is unbelievably quiet, which my patients love. I would definitely find it difficult to work without it!"
Alina Sheikh BDS, Principal Dentist at Enhance Dental Care, East Kilbride, www.enhancedentalek.co.uk

S·Max M

- Stainless steel body
- 4 hole spray
- 2 head sizes - standard and miniature
- DIY cartridge change

Buy 3
& get
extra 1
FREE*
(optic)

Buy 3
for
£1,095*
(non-optic)



<div><div>Ti-Grade</div><div>Optic</div></div>		S-Max Advanced Turbines				
		Fitting	Model	Order Code	RRP	Warranty
		NSK	M900L	P1254001	£729	2 years [†]
		KaVo®	M900KL	P1258001		
		Sirona®	M900SL	P1262001		
		W&H®	M900WL*	P1266001		
		Bien-Air®	M900BL*	P1264001		
Non-optic	NSK	M900	P1256001	£679		
	KaVo®	M900K	P1259001			
	W&H®	M900W	P1287001			

* Compatible Turbines with Integral LED

Pana·Max2

Buy 1
get 1
FREE*



- Save time and money - DIY cartridge change
- Enhanced power and endurance compared to previous models††
- Slimline, scratch-resistant stainless steel body
- Reinforced head ensures durability

††30% more power

Pana-Max2 features a reinforced head with NSK bumper technology for added longevity and protection of the cartridge.



Non-optic	Pana-Max2 Budget Turbines					
	Pana-Max2	Fitting	Head size	Order Code	RRP	Warranty
	Direct connection + single hole spray	B2 (Borden)	Ø 11.2 x H 13.4mm	P1155001	£372	12 months
		M4 (Midwest)		P1154001		
	Coupling connection + single hole spray	PTL (NSK)		P1211001	£ 442	
		KV (KaVo®)		P1212001		

For the complete range of contra-angles visit www.NSK-uk.com

*Terms and conditions apply. Offers only available through participating dealers.
**One per instrument



Handpiece oil – Pana Spray Plus
Z182100 or Y900630 RRP £30.70

Buy 1 get 1 FREE*

For further info on handpiece care visit www.nskcare.co.uk

VIVA ace

Ideal as a
surgery
back up

Monthly
rental
Only

£433*

Fully maintained
for 3 years

**FREE S-Max
M95L, M25L
& M65
handpieces***

Height
28.4 cm
Depth
29.8 cm
Width **46.2 cm**



Mobile Dentistry System

- Portable (only 8.6kg) dental treatment centre
- Lightweight micromotor with additional endo functions
- Multipurpose ultrasonic scaler
- Quick setup and storage
- 4-level adjustable height and angle

Y1003773
VIVA ace Complete Set
RRP £11,346

*Terms and conditions apply. Offers only available through participating dealers.

Air Polishers

Prophy-Mate neo

Air Polisher



Buy Prophy-Mate neo and get FREE second handpiece worth £378*

Supragingival Air Polisher

- Compact and lightweight
- Includes 60° and 80° nozzles
- 360° rotation
- Packed in a handy case
- Easily dismantled for autoclaving
- Includes 6 sachets of FLASH pearl® that quickly removes stains and plaque

Prophy-Mate neo RRP £1018.80

Y135029 for NSK PTL coupling

Y135030 for KaVo® MULTiflex® LUX coupling

Y135031 for Sirona® quick coupling

Y135032 for W&H® Roto Quick® coupling

Y135033 for Bien-Air® Unifix® coupling

Prophy-Mate neo nozzle

T797012 60° nozzle RRP £109

T797013 80° nozzle RRP £109

Buy 1 get 1 FREE*

Find out more at www.myNSKoralhygiene.co.uk



Buy 2 packs get a Prophy-Mate neo nozzle FREE* SAVE £106

4 x 300g bottles
Y900698
RRP £125

100 x 15g sachets
Y900693
RRP £171

Perio-Mate

Biofilm Eraser



SAVE £452 Buy both Perio-Mate & Prophy-Mate neo for £1,675*

Subgingival Biofilm Eraser

- Avoids injury to delicate periodontal pockets
- Smooth 360° rotation
- Flexible plastic nozzle tip for better access to difficult to reach areas
- Ultra-fine water spray and Perio-Mate powder eradicate biofilm in just 20 seconds
- Packed in a handy case

Perio-Mate RRP £1,161

Y1002654 for NSK coupling

Y1002655 for KaVo® MULTiflex® LUX coupling

Y1002656 for Sirona® quick coupling

Y1002657 for W&H® Roto Quick® coupling

Y1002658 for Bien-Air® Unifix® coupling

Find out more at www.myNSKoralhygiene.co.uk



Perio-Mate Powder
RRP £87
Glycine based powder,
2x 160g bottles
Y900938

Perio-Mate Tip Nozzles
RRP £90
Pack of 40
Y1002741

"I use the Perio-Mate routinely as it is ergonomically designed, lightweight and offers very good balance in my hand. Additionally I like the disposable tips as they reduce maintenance. Perio-Mate, with its ultra-fine water spray and glycine-based powder, is part of my standard implant maintenance protocol for mucositis and early peri-implantitis as it delivers excellent results and patients find the procedure surprisingly painless and the taste rather pleasant."

Dr. Chong Lim. BDS, MSc, MSc, Implantologist & Specialist Periodontist, Perio & Implant Dental Clinic Referral, Twickenham, London, www.perionimplant.com

Prophy Handpieces

FX57m

FX65m


AR-EC



NAC-EC
THE ORIGINAL



NAC-EC Triple pack offer £180*



Buy 1 FX get 1 FREE*

	Model	Max speed	Brushes / Cups	Code	RRP
FX Series	FX57m	10,000rpm	Screw in	C1054001	£270
	FX65m	40,000rpm	Disposable Prophy heads	H1014001	£239
EX Series	AR-EC(s)	5,000rpm	Screw in	Y110164	£165
	AR-EC(k)	5,000rpm	Snap on cups	Y110165	£165
	NAC-EC	30,000rpm	Latch grip	Y110148	£141

*Terms and conditions apply. Offers only available through participating dealers.

Varios Combi Pro

Two oral care systems in one!



FREE
Perio Set worth
£1,161.30 + VAT*‡



OR

Half price
for NSK
Prophy-Mate
or Surgic Pro
users*‡

Varios Combi Pro RRP £4,511
Combined Oral Hygiene Unit
Y1002843

Two systems in one - NSK has successfully combined ultrasonic treatment and air polishing into a single professional treatment unit that streamlines procedures and saves time and space.

Note: Connection to external 6mm air & water line sockets are needed to run this unit

- Automatically sets air polisher functions for supragingival or subgingival applications
- Independent unit and chamber for easy switching between applications
- Multi-foot control delivers all key operations
- Optional perio system also available
- Wide range of applications
- Gentle on soft tissue
- Wide choice of tips
- A clear functional display
- Auto-cleaning function
- Easy maintenance

NSK



IKIGAI

ORAL HYGIENE PROGRAMME

Working with



FIND YOUR IKIGAI

The NSK Ikigai Oral Hygiene Programme offers a new way of working to deliver exceptional oral care coupled with a deeper level of job satisfaction - enhancing your practice and your patients' wellbeing.

- More confidence, energy and motivation to develop your career and oral care practice
- A deeper understanding of current oral hygiene principles and how to maximise earnings
- A greater sense of fulfilment by being in control of your daily practice
- An understanding of how to reduce stress and frustration
- Clarity about how to build more effective working relationships

6 HOURS CPD
+ a goody bag
worth £75!



生き甲斐

FIND YOUR IKIGAI - BOOK AT
www.nsk-ikigai.co.uk

*Terms and conditions apply. Offers only available through participating dealers. ‡ The 2 offers above cannot be combined.

Varios 970

iPiezo engine





Call 0800 6341909
for details



Buy Varios 970 /
Varios 970 LUX and
get 2 additional
VA2 handpieces
and sterilisation
cassettes FREE*

SAVE up to
£1,002.50
+ VAT*


Ultrasonic Scaler

Ultrasonic scalers not only offer patients a more comfortable experience as less water is needed, but offer the user a wide range of 70 tips, not only for scaling/perio/implant maintenance, but for endo, extractions, restorative & MI applications.

Varios 970 LUX RRP £2,627 Optic Complete set Y1004161	Varios 970 RRP £2,260 Non-optic Complete set Y1004160	VA2-LUX-HP RRP £482 Varios optic handpiece E351050	<div>TRIPLE PACK OFFER*</div> <div>OPTIC £999 NON-OPTIC £585</div>	VA2-HP RRP £269 Varios non-optic handpiece E350050
---	---	--	--	--

Buy 1 get 1
free on all
scaler tips*

NSK Ultrasonic Scaler Tips



Find out more and see the full range of Varios tips at www.myNSKoralhygiene.co.uk



Perio/Hands-on courses proudly
supported by NSK
www.periowithdeepak.com



*Terms and conditions apply. Offers only available through participating dealers.

Surgic Pro

Surgical Implant Unit



2 year warranty

FREE
2 additional
20:1
contra-angles
with each
model*

Also FREE
with model C:
Osseo 100*
see below

Buy extra
SurgicPro
motor/cord get
20:1 implant
contra-angle
FREE*

Surgical Implant Unit

Surgic Pro

Choice of 3 models. Each model comes with a control unit, motor/cord, 20:1 implant contra-angle and accessories.

MODEL A Surgic Pro non-optic RRP £4,062
with S-Max SG20
Y1003587

MODEL B Surgic Pro LED RRP £5,441
with Ti-Max X-SG20L
Y1003586

MODEL C Surgic Pro+ LED RRP £6,054
with Ti-Max X-DSG20L
Control unit has USB port for recording all operations
Y1003585

Accessories

SG70M RRP £1,679.90
Surgic Pro non-optic motor/cord. Comes with 'Model A'
E1025001

SGL70M RRP £2,416
Surgic Pro LED motor/cord. Comes with 'Model B' & 'C'
E1023001

Irrigation Tubing RRP £57
Pack of 5. Comes with all options
Y900079

"I have found NSK's Surgic Pro very easy to use for surgical and implant procedures. It offers powerful and precise performance, coupled with intuitive operation and has revolutionised my implant surgery."

Dr Joe Bhat BDS FDS RCS MCLinDent MRD RCSEd, Specialist in Prosthodontics and Oral Surgery

Already a user? Contact us
about our peace of mind
service contract options for
NSK surgical equipment.

Introducing Osseo 100

Osseointegration Monitoring Device


Measure implant stability and osseointegration to enhance decisions about when to load the implant. Especially important when using protocols with shorter treatment time and treating higher risk patients.

Osseo100 RRP £2,075
Code Y1004175

30% off
Only
£1,450*

Speak to your local NSK representative today to find out more

3-step procedure



Implant Contra-angles 20:1

- Max speed: 2,000min⁻¹
- Warranty: S-Max
18 months, Ti-Max 2 years

TI-MAX



S-MAX



For a full range of our implant/
surgical handpieces, see
www.nsk-uk.com

Optic Ti-Max X-DSG20L
RRP £1,385
Dismantleable.
Comes with 'Model C'
C1068001

Optic Ti-Max X-SG20L
RP £1,155
Comes with 'Model B'
C1003001

Non-optic S-Max SG20
RRP £975
Comes with 'Model A'
C1010001

Surgical Straights 1:1

- Used with standard HP burs
- 2 year warranty



TRIPLE
PACK
OFFER*
NON-OPTIC £1185
OPTIC £1680

Non-optic
Ti-Max X-SG65 RRP £577
H1038001

Optic
Ti-Max X-SG65L RRP £846
H1009001

*Terms and conditions apply. Offers only available through participating dealers.

VarioSurg3



VarioSurg3
RRP £7,330
with LED handpiece
Y1002726

Ultrasonic Surgery System with LED

- Piezo ultrasonic bone surgery system
- Intuitive Piezo feedback
- Large, clear LCD screen
- Wide range of tips ensures faster and more precise treatments
- Lightweight, slimline VarioSurg3 handpiece, with 50% more power than the VarioSurg2

+ FREE 2nd
handpiece/cord or
12 additional tips

OR

50% off for existing
Surgic Pro/XT
users* ‡

"I've been using NSK Surgical equipment for the past few years and have found it to be exceptionally durable, well designed and predictable to use. I am hugely appreciative of and impressed by NSK's customer service."

Dr. Hugh Byrne BDS, Dalkey,
County Dublin.
www.dalkeyclinic.com

Surgic Pro+ × VarioSurg3



Surgical Link Stand

Link together your Surgic Pro and VarioSurg3

Control both units with
one foot pedal. Choose
a Link Stand or store
everything on an iCart
where you can adjust
the height and move
around on its wheels.

SG Link Stand Set RRP £465
Y1002729

iCart Duo RRP £653
S9090



iCart Duo

Buy any
Surgic Pro
and get
VarioSurg3
HALF
PRICE*

Buy
Endo-Mate
DT and get
iPex II FREE
worth £898*

Endodontics

- Large LCD display
- 5 programmes for different file systems
- Supports most major brands of Ni-Ti files
- Torque control facility ranging from 0.1 to 4.0 Ncm depending on head:gear ratio
- Auto-reverse and alarm function sounds an audible alert when preset torque level has been reached

Endo-Mate TC2 RRP £1,549.80 Y1001027

ENDO-MATE TC2



ENDO-MATE DT



- Large LCD display
- Store up to 9 speed and torque settings
- Use with Ni-Ti files from most major suppliers
- Torque control from 0.1 to 6.5 Ncm depending on head:gear ratio
- Auto-reverse on, auto-stop, auto-reverse off

Endo-Mate DT
RRP £1,465
Y141121

iPexII APEX LOCATOR



- Accurate measurement of any root canal using SmartLogic technology
- Operates in all conditions, including dry, wet and bloody canals

iPexII RRP £921 Y1002208

*Terms and conditions apply. Offers only available through participating dealers. ‡ The 2 offers above cannot be combined.

iCare

Handpiece oiling & maintenance unit

**Choose
50% off or
FREE handpieces
worth £1,937***

2x S-Max M900L turbines,
1x FX25m contra-angle
and 1x FX65m straight
handpiece*

**Handpiece Oiling Unit**

- Standardised oiling reduces handpiece repair bills
- Works the gears of speed increasing handpieces
- Oils up to 4 instruments (10 seconds each)
- Size: H350 x W280 x D240mm

Note: Requires fitted 6mm air line at 5-6 bar

iCare C2 RRP £2,330

2 high speed, 2 slow speed ports
Y1002787

iCare C3 RRP £2,448

1 high speed, 3 slow speed ports
Y1002788

iCare+

Handpiece decontamination unit

**SPECIAL BUY
£3,999 or
£5,618 including
FREE handpieces
worth £4,307***

2x Ti-Max Z900L turbines,
1x Z25L contra-angle
and 1x X65L straight
handpiece*

**Handpiece Decontamination Unit**

- Cleans, disinfects and lubricates handpieces (like a washer disinfectant for handpieces)
- For up to 4 instruments
- Size: H355 x W400 x D405mm

Note: Requires fitted 6mm air line at 5-6 bar

iCare+ C2 RRP £5,618

2 high speed, 2 slow speed ports
S102001

iCare+ C3 RRP £5,618

1 high speed, 3 slow speed ports
S103001

iClave plus

B-type vacuum autoclave

**Only
£3,995
saving £1,623***



- 20% more sterilisation space than a conventional stainless steel chamber of the same size
- Copper chamber minimises internal temperature fluctuations
- Super fast cycle as quick as 18 min, including drying
- 2 year warranty
- Unit size: H428 x W445 x D532mm, chamber 18L

MODEL iClave plus 230V RRP £5,618
REF **DXY1003077**

iMax S

S-type autoclave

**Only
£2,395
saving £423***



- 28 minutes for an S cycle
- 20 minutes for the Rapid S class cycle. (times inc drying)
- Rapid cycles for unwrapped or solid loads

- 18 litres – 20% more capacity†
- Copper chamber provides superior thermal conductivity
- Low energy consumption
- USB data logger
- Efficient drying phase, gentle for handpieces
- 1 year warranty
- Separate reservoirs – doesn't re-use water
- Easy to use control panel

MODEL iMax S RRP £2,818
REF **DX91803S0020**

†Than conventional autoclaves.

*Terms and conditions apply. Offers only available through participating dealers.