

NSK

Pana Max 2

NSK Pana Max 2
high speed handpiece
for connection directly
on to your existing tubing.
Single Spray.
Non Optic
Midwest 4 Hole /
Borden 2 Hole
12 Months Warranty



www.ukdentalsupplies.com/nsk

Pana Max 2 PTL/KV

NSK Pana Max 2
high speed handpiece
for connection on to
NSK PTL or Kavo Multiflex
couplings.
Single Spray
Non Optic
12 Months Warranty



NEW & IMPROVED MORE POWER S-Max M900

NSK S Max M900 turbine
non optic / M900L Optic.
Torque / Standard head.
4 Port Spray
Available to fit NSK,
Kavo, W&H couplings.
26 watt power output.
2 Years Warranty



M900 TRIPLE PACK
NSK OFFER: £1,095.00

BDSI PRICE: £876.00

M900L RRP: £729.00

BDSI PRICE: £524.00

BUY 3 GET 1 FREE

KAVO Kerr

www.ukdentalsupplies.com/kavo

SMARTtorque Turbine

Kavo Smarttorque S619L Optic
S619C Non Optic
16 watt power output
Non-Optic
RRP: £542.00

BDSI PRICE: £315.00

Optic
RRP: £677.00

BDSI PRICE: £365.00



EXPERTtorque Turbine

Kavo Experttorque E680L / E680C
18 watt power output
Connection to Kavo couplings
4 Port Spray
18 Months Warranty

Non-Optic
RRP: £980.00

BDSI PRICE: £479.00

Optic
RRP: £1,047.00

BDSI PRICE: £515.00

QUAD PACK PRICE £1,545.00



MASTERtorque Turbine

One of the quietest high-speed
turbines on the market.
Output of up to 23 watt
Direct Stop Technology -
1 second stop.

RRP £1,430.00

BDSI PRICE: £749.00

BUY 4 GET 1 FREE




**BUY PRODUCTS
AND SERVICES
FROM SUPPLIERS
YOU CAN TRUST**

BDIA Member companies
Suppliers you can trust



BRITISH DENTAL
INDUSTRY
ASSOCIATION



MK-dent

High quality dental turbines manufactured in Germany
www.ukdentalsupplies.com/mkdent



ECO LINE

MK-dent HE21KL
24 watt power output
Compatible with NSK, Kavo,
W&H, Sirona Couplings.
18 Month Warranty

NON OPTIC: £249.00
OPTIC: £339.00



CLASSIC LINE

MK-dent HC20KL
Classic Line Power Edition
28 watt power output.
Compatible with Kavo
couplings only.
18 Month Warranty

NON OPTIC: £429.00
OPTIC: £515.00



PRIME LINE

MK-dent HP21L
Vibration free running.
Sound level less than 55dB.
25 watt power output.
Compatible with NSK,
Kavo, W&H, Sirona Couplings.
24 Month Warranty

OPTIC: £449.00

15% DISCOUNT WHEN YOU BUY 3 MK DENT TURBINES

glows™
Dental Equipment. Taiwan

Super-Max Turbine

Glows Super-Max Turbine
4-hole midwest connection direct to tubing
4 port water spray
Easy DIY cartridge replacement
18w power output
Made in Taiwan
12 month warranty

£149.00

£399.00 TRIPLE PACK



Super-Max Optic Turbine

Glows Super-Max Optic turbine - NSK connection
Glass fibre optic rod
4 port water spray
Easy DIY cartridge replacement
18W power output
Made in Taiwan
12 month warranty

£250.00

£695.00 TRIPLE PACK



ND NAKAMURA DENTAL

www.ukdentalsupplies.com/nakamura

High Quality Dental Handpieces Manufactured in Tokyo, Japan.



ME-20M Air Motor 20,000RPM
24 Months Warranty

£150.00



ECN-20L Latch Grip Contra Angle NSK Compatible
12 Months Warranty

£65.99

BUY 3 @ £49.00 each



ES-30A Straight 1:1
12 Months Warranty

£109.00

BUY 3 @ £98.00 each

KAVO Kerr

Kavo SMARTmatic S20

Contra Angle 1:1
Max speed 40,000rpm.
For CA burs

RRP £208.00

BDSI PRICE: £127.00



Kavo SMARTmatic Endo S81

Contra Angle 8:1
Miniature head for
Endodontic procedures

RRP: £300.00

BDSI PRICE: £182.00



Kavo SMARTmatic S10

Straight Handpiece 1:1
Max speed 40,000rpm
For HP burs

RRP: £200.00

BDSI PRICE: £124.00



NSK

FX-25M

NSK FX-25M 1:1 Direct drive push button chuck TD. Longer neck and smaller head. Max speed 40,00rpm.

RRP: £233.00

BDSI PRICE: £225.00

BUY 1 GET 1 FREE



FX-15M

NSK FX-15M 4:1 reduction push button chuck TD. Longer neck and smaller head. Max speed 10,000rpm

RRP: £269.00

BDSI PRICE: £255.00

BUY 1 GET 1 FREE



NAC-EC

NSK NAC-EC Non-optic latch type contra angle handpiece for use with CA/RA burs. Max speed 25,000rpm

RRP £137.00

BDSI PRICE: £79.50

TRIPLE PACK £180.00



MK-dent

AM1004

Airmotor with Internal Water 20,000 rpm 4-hole connection Spray regulation system 2 year warranty

£249.00



AM1116

Airmotor with internal water and LED optics 20,000rpm 6-hole connection Spray regulation system 2 year warranty

£339.00



LE11

MK Dent Contra Eco Line 1:1 LE11 - Non optic LE11L - Optic internal water 40,000rpm 2 year warranty

Non Optic £249.00

Optic £299.00



Bien Air[®]

Dental

OPTIMA PORTABLE MICROMOTOR SYSTEM

Equipped with the MCX-LED Brushless Micromotor with 3 year warranty. 2 year warranty on main Optima unit. Forward / reverse operation for increased flexibility and versatility SmartLogic control for adjustable speed. Internal Water. LED light for optimal illumination. Brushless technology for extended service life Speed range: 1,000-40,000 rpm Maximum torque: 2.5 Ncm



CALL TODAY FOR A FREE TRIAL!



Bien Air RRP: £1,541.50

BDSI PRICE: £985.00

Micro series handpieces shorter and lighter with 2 year warranty.

CA 1:1L

Dualook 2 optical glass conductors Speed 600 - 40,000 rpm.

RRP £733.59

BDSI PRICE: £390.00



CA 1:5L

Increasing speed can be used for caries removal. Dualook 2 optical glass conductors. Speed 600 - 200,000rpm.

RRP £902.48

BDSI PRICE: £480.00



MC3

Electric Micromotor Non-Optic

RRP: £736.00

BDSI PRICE: £449.00

Optic

RRP: £920.00

BDSI PRICE: £560.00



More available on www.ukdentalsupplies.com
call our fully trained staff for advice - 01226380400

NEW Table-Top Implant AND Oral Surgery Unit from Bien Air
Carry out oral surgery procedures, particularly wisdom tooth extractions, when combined with the PM 1:2 straight handpiece
Perform oral and periodontal procedures when combined with the new CA 1:2.5 high speed handpiece.
2 year warranty on unit - 3 year warranty on MXi PLUS motor
RRP: £3,017.97



BDSI PRICE: £2,095.00

Bien Air PM 1:2

1:2 speed increasing straight handpiece for implantology
Ideal for wisdom tooth extraction
Speed range: 200-80,000 rpm
External irrigation on detachable sleeve
Bur locking ring system
1 year warranty



RRP: £1,150.91

BDSI PRICE: £575.00

Bien Air CA 20:1 L

Features small head for implant work. Irrigation line routed inside handpiece
Maximum speed: 2,000 rpm
Head diameter: 9.8mm
Head size: 13.5mm
Weight: 80g, Noise level: 57 dBA
2 year warranty



RRP: £1,229.00

BDSI PRICE: £845.00

Bien Air CA1:2.5L

For a wide range of applications;
Wisdom tooth extraction,
Hemisection Apicectomy,
Crown Lengthening.
Irrigation Line routed inside handpiece
1 year warranty



RRP: £1,350.00

BDSI PRICE: £795.00

NSK SURGIC PRO

The NSK Surgic Pro is compact, lightweight and powerful
SGL70M micromotor is 16.2mm shorter and weighs 42g less compared to the conventional optic surgical micromotor.
The unit can memorise eight different implant systems and a 64 programmes.
The programmable parameters are gear ratio, speed, rotation direction, torque limit, coolant solution volume and illumination intensity.



Model A; Non Optic Surgic Pro with SG20 Handpiece RRP £4,062.00 **BDSI PRICE: £2,924.00**

Model B; Optic Surgic Pro with TI Max X-SG20L RRP £5,441.00 **BDSI PRICE: £3,917.00**

SG20

20:1 reduction implant surgical handpiece. Non-Optic
Max Speed 2,000 min-1
Max Torque 80 Ncm



Non Optic
RRP £975.00

BDSI PRICE: £701.00

Ti-Max X-SG65L

Straight 1:1 drive surgical handpiece.
Accepts HP burs (Ø2.35).
Also accepts CA burs (Ø2.35).
External cooling.
Maximum speed 40,000 rpm



Non Optic RRP £577.00

BDSI PRICE £415.00

Optic RRP £845.60

BDSI PRICE £794.00

Ti Max Z-SG45L

1:3 Speed Increasing Surgical Handpiece with external irrigation
45° angled head doesn't interfere with anterior teeth, giving better access to the treatment area
3 year warranty.



Non Optic RRP £879.40

BDSI PRICE £633.00

Optic RRP £1,105.90

BDSI PRICE £796.00

ISD900 PROSTHETIC SCREWDRIVER

Cordless Implant Prosthetic Screwdriver for insertion of screws and implants.
ISD900 can be used with any dental implant system to accurately tighten implants to the precise torque value as per the manufacturers specification.
Torque 10-40Ncm in 1 or 5 Ncm increments

RRP £1,436.00

BDSI PRICE: £1,033.00



NSK

PROPHY-MATE NEO



Air powered tooth polishing system.
Includes 6x Flash Pearl Sachets.
To fit NSK, Kavo, Bien Air, W&H, Sirona couplings.

RRP: £994.00 **BDSI PRICE: £733.00**

PERIO-MATE



NSK Perio-Mate biofilm eraser
Avoids injury to delicate periodontal pockets.
Smooth 360° rotation
Flexible plastic nozzle tip for better access to difficult to reach areas
Ultra-fine water spray and Perio-Mate powder eradicate biofilm in just 20 seconds
To fit NSK, Kavo, Bien Air, W&H, Sirona Couplings.

RRP: £1,161.00 **BDSI PRICE £836.00**

VARIOS 970



Available in Optic / Non Optic
Increased Patient comfort
wide range of tips available
for scaling, endo, extractions etc

BDSI PRICE NON OPTIC: £1,627.00

BDSI PRICE OPTIC £1,891.00

PLUS GET 2 x VARIOS HANDPIECE FREE

FLASH PEARL POWDER

100 x 15g Sachets
RRP: £171.30

BDSI PRICE: £123.34



PERIO-MATE POWDER

Glycine based powder
2 160g bottles
RRP: £90.00

BDSI PRICE: £69.68



**BUY PROPHY-MATE NEO &
PERIO-MATE TOGETHER FOR
£1,255.00**

VARIOS SCALER HANDPIECES



VA2-HP Varios Non- Optic Handpiece

BDSI PRICE £215.00

TRIPLE PACK PRICE £585.00

VA2-LUX-HP Varios Optic Handpiece

BDSI PRICE £385.00

TRIPLE PACK PRICE £999.00

ACTEON

NEWTRON BOOSTER



Connect direct to your water supply
Uses Newtron handpiece (F12281) - non LED
The tip colour coding system
associates each tip with each of the
4 power ranges. 5 year main unit warranty

RRP £1,081.50 **BDSI PRICE £815.00**

NEWTRON P5 XS

With Old Newtron Handpiece



With indepent water reservoir
powered by a peristaltic pump.
Uses Newtron Handpiece(F12281)- non LED
The tip colour coding system associates each tip with
each of the 4 power ranges. 5 year main unit warranty

RRP £1,705.00 **BDSI PRICE £1,279.00**



Bonart Dental

BONART E1 ELECTROSURGERY

10 power level settings. 3 operating modes:

Cut, Coag 1 & Coag 2.

Includes 7 piece electrode set and indifference plate

Monopolar unit. 12 months warranty

Very popular, durable unit

£495.00



Bonart ART M3II

25K / 30K

Cavitron compatible Scaler

BEST SELLER

£349.00 + 1 Insert FREE



Bonart Scaling Inserts

compatible with all other brands

25/30k. Universal, Slimline,

Beavertail & external spray

£46.00

BUY 5 GET 1 FREE



Bonart ART M4Pro

25K / 30K

Water bottle supply

£495.00

+2 Inserts FREE

DTE



DTE D7

Ultrasonic Piezo Scaler -

Complete with 6 tips and

Satelec compatible handpiece

£270.00

£295.00 (LED model)



DTE D5

Ultrasonic Piezo Scaler -

Complete with 6 tips and

Satelec compatible handpiece

£164.00

£185.00 (LED model)



Woodpecker UDS-E

Ultrasonic Piezo Scaler -


Complete with 5 tips and

EMS Compatible handpiece

£325.00

£345.00 (LED model)

Replacement Scaler Handpieces. DTE / Woodpecker. Satelec / EMS Compatible



£89.95 BUY 2 GET 1 FREE



Replacement tips

£12.50

BUY 5 GET 1 FREE



Woodpecker X-Cure

1 Second Cure. High Power
LED Curing Light.
Output wavelength 385nm-515nm

£295.00



Woodpecker LED-C

Portable LED curing light with
2 interchangeable LI-Ion batteries.
Output wavelength 420nm-480nm

£105.00



Saeyang K38

Portable Micromotor kit
Control unit, electric motor,
Straight Handpiece

£299.00



Automatic Crown & Bridge Remover

Anthogyr Automatic Crown and Bridge Remover
Fastest on the market (up to 20 strokes/second)
the most efficient and least stressful.

£465.00

- 1 x Long angulated hook (1428)
- 1 x Short angulated hook (1427)
- 1 x Flat hook (1429) non slip grooves
- 1 x Wire tip ATD 12lg 50mm
- 1 x Wire tip ATD 12lg 60mm
- 1 x Ceram wire-tip ATD 22lg 95mm
- 1 x spare wire lg 50mm
- 1 x spare wire lg 60mm
- 1 x spare wire lg 150mm



Linea-tac TAC 400M

Digital Amalgamator
Made in Italy

£195.00



Hygea D3ntal

Ultrasonic Bath Compliant to
HTM01-05 and ISO standards

£1,675.00



MX -300

Alginate Mixing Machine
produce perfectly mixed alginate
at the touch of a button.

£675.00

Acteon PSPIx2 Imaging Plate Scanner



Create digital dental x-rays in seconds with the PSPIx 2 X-ray scanner.
Touch screen operation with simple to use interface. Can work on up
to 10 workstations- perfect for multi surgery dental practices.
Fully automatic processing - insert the phosphor plate into the PSPIX2
and no further input is needed. the image will be processed in a few seconds.
ACTEON RRP: £6,484.00. **BDSI PRICE: £4,573.00**

PSPIx 2 Phosphor Plates



Pack of 2

- Size 0 - 22 x 31 mm
- Size 1 - 24 x 40 mm
- Size 2 - 31 x 41 mm
- Size 3 - 27 x 54 mm

£60.00



GP CUTTER

Gutta Percha Cutter with three cutting tips. Simple 1-button operation heats the tip to 500 °c in 1 second

£90.00



DPEX III GOLD

Apex Locator features LCD screen. 97.71% accuracy.

£149.00



Endo A Class LED

Cordless endodontic motor with LED. Reciprocating. 9 Memory Settings, Li-Ion Battery. Speed 60rpm - 900rpm

£475.00

ENDODONTIC FILES

FULL ENDO RANGE AVAILABLE VISIT www.ukdentalsupplies.com



Multi Taper Files - Pro-taper System

SX | S1 | F1 | F2 | F3 - Pack of 6

£16.99



Greater Taper Files - Rotary System

.4 Taper - 15 | 20 | 25 | 30
35 | 40 - Pack of 6

.6 Taper - 15 | 20 | 25 | 30
35 | 40 - Pack of 6

£11.99



Wave + Gold - One File System

Small | Primary | Medium | Large - Pack of 6

£19.99

DENTAL STOOLS MANUFACTURED IN THE UK



SADDLE STOOL

Dentist Saddle Stool with tilt mechanism
Ergonomically designed to retain correct posture

£320.00



SURGEON STOOL

Surgeon Stool with deep contoured back rest and seat for extra comfort

£310.00



ASSISTANT STOOL

Back Support serves as armrest which can be positioned by the factory on the left, right or in the middle

£275.00

Visit the website to view the full range. Upholstered in Naugahyde Spirit Millenium vinyl fabric.

All prices exclude VAT. Offers valid until 30/04/2020. Subject to availability, e&o excepted. Offers cannot be combined

www.ukdentalsupplies.com | 01226 380400 | sales@bdsi-ltd.co.uk