



KAVO

Dental Excellence

**Handpieces &
Small Equipment**

Special Offers

BDSI

www.ukdentalsupplies.com

60
YEARS OF KAVO UK
Est 1962

Hello!

We're delighted to present the highlights of our current Handpiece and Small Equipment promotions. All of our most popular handpiece models are available at a special promotional price and some have added value offers too.



KaVo Turbines



Entry Level

SMARTtorque

12 MONTH WARRANTY

TRIPLE PACK

£833.00 Triple Pack Price



Expert Level

EXPERTtorque

18 MONTH WARRANTY

DUO PACK

£959.00 Duo Pack Price



Master Level

MASTERtorque

24 MONTH WARRANTY

Direct Stop Technology™

PlasmaTec Coating™

3+1 available

£845.00 each

KaVo ELECTROmatic

The KaVo **ELECTROmatic** converts your treatment unit's air power into electric power to facilitate the use of speed increasing handpieces.

We have several versions: with short cables, long cables, dual motors and endo function – and all of them are **simple to install**.

We've got a **special offer** on the **ELECTROmatic Plus** with **three speed increasing handpieces** of your choice at a very special price.



ELECTROmatic Plus, with KL703 motor and 3 x Speed Increasing handpieces

Package Prices from £2,995.00

Speed Increasing Handpieces

KaVo's speed increasing handpieces are designed to perform a wide range of clinical functions to the highest levels of performance.



KaVo MASTERmatic M25L 1:5
24 MONTH WARRANTY

Available as part of a
buy 3 get 1 free bundle deal.

from £854.00 each



KaVo EXPERTmatic E25L 1:5
18 MONTH WARRANTY

Available as part of a **DUO Pack**.

£1,195.00 DUO Pack Price

KaVo RONDOflex Plus 360

Our superior air abrasion unit with additional water supply and reduced powder load for a wide range of minimally-invasive applications. With a **KaVo MULTIflex** coupling, the **RONDOflex** is fully transportable and can be used in any practice.



RONDOflex Plus: **£1,100.00**



Scan the **QR code**
for more details and
case studies.

KaVo PROPHYflex 4 & PERIO Kit

The KaVo PROPHYflex 4 air polishing system offers superior cleaning - removing stains, softening calculus, cleaning pockets and leaving patients' teeth feeling and looking cleaner.

The PROPHYflex 4 has three power settings to cater for even the most sensitive patients, an optional Power Tip, and a special adapter Kit which will allow you to use it for subgingival cleaning. This all-rounder comes in an added value pack with our new generation PROPHYpearls for enhanced cleaning that is both gentle and highly effective.

PROPHYflex 4: from £860.00



Scan the **QR code**
for more details and
case studies.



QUATTROcare & SPRAY

The KaVo QUATTROcare cleans and oils up to four handpieces at once in under a minute. The automated, metered dosing saves you time and money on cleaning, oiling and maintenance.

QUATTROcare: **£1,365.00**


Plus 6 x cans of KaVo QUATTROcare Spray FOC!



SONICflex & Tips

The KaVo SONICflex is the next generation air scaler designed for ease of use and highly efficient scaling. Available with the widest range of tips for unparalleled versatility and convenience.

SONICflex 2003L **£985.00**



Surgical Handpieces

We have a full range of straight and contra-angle surgical handpieces for cutting-edge performance.

from **£755.00** or get **one FREE** with the **MASTERSurg** or **EXPERTSurg**.



MASTERSurg & EXPERTSurg

The KaVo surgical units offer programmable settings, a colour screen and an auto-calibration feature to test the motor and the drill quickly and easily before each use.

Popular in hospital settings as well as high street clinics, choose between the MASTERSurg - with a touch screen and wireless foot control - or the EXPERTSurg. Each one comes with a FREE KaVo surgical handpiece!

MASTERSurg / EXPERTSurg: from **£2,869.00**

