

Autumn/Winter Offers



NEW

VALID UNTIL
31/12/2023

Power that drives performance

**Unparalleled power,
performance and protection.**

The new Ti-Max Z990L and Z890L Turbines fuse the very best of NSK technology, making treatment faster, safer, more comfortable and stress-free - for clinicians *and* patients.

2 + 1 FREE

+ a further **FREE Z micro**
with ultra mini head.



Power that drives
performance
from 30- 44W*




*Ti-Max Z990L with air pressure at 3 to 4.2 bar

Design that delivers
performance
100° angle head



Built to
deliver
performance
**lightweight
and durable**



**SCAN FOR MORE
INFORMATION**

Innovation
that defines
performance
**New Quick-Stop
function**




NEW Ti-Max Z Series High Performance Turbines



Unprecedented power - The new "DYNAMIC POWER SYSTEM" delivers between 30-44W** of power, transforming time-consuming procedures involving even the hardest of materials, including zirconia.

*Ti-Max Z990L with air pressure at 3 to 4.2 bar.

100° angle head - Designed to work with the curves of the oral cavity, the new head angle allows easier access to the molar area with better visibility, resulting in more efficient and effective treatment.



New Quick-Stop function - High-speed rotation can be halted in an instant, reducing the risk of suckback, providing reliable protection against cross-infection and ensuring safer treatment.

Lightweight and durable - The titanium bodies of these new NSK turbines are even lighter than their predecessors and, combined with DURAGRIP, provide a supremely comfortable grip to reduce hand fatigue, resulting in more efficient and effective treatment.



Ti-Max High performance turbines



2 + 1 FREE* + 1 FREE Z micro

NEW

• Lightweight external titanium body • DYNAMIC POWER SYSTEM • 3 year warranty



- Optic
- Power 30-44w**
- Std head: Ø 12.1 x 12.7 mm

MODEL Z990L
CODE PA23730001
RRP £1,134

Also available in fittings for
KaVo®, Dentsply Sirona®, W&H®, Bien-Air®

**At 3 - 4.2 bar



- Optic
- Power 28-41w
- Mini head: Ø 10.6 x 10.8 mm

MODEL Z890L
CODE PA23720001
RRP £1,134

Also available in fittings for
KaVo®, Dentsply Sirona®, W&H®, Bien-Air®

For full info on
models/order codes

Ti-Max



1 + 1 FREE*

• Lightweight external titanium body • Clean head system • 4 hole spray • 3 year warranty



- Optic
- Power 26w
- Std head: Ø 12.5 x 13.1 mm

MODEL Z900L
CODE P1111001
RRP £1,134

Also available in fittings for
KaVo®, Dentsply Sirona®,
W&H®, Bien-Air®



- Optic
- Power 23w
- Mini head: Ø 10.8 x 12.1 mm

MODEL Z800L
CODE P1110001
RRP £1,134

Also available in fittings for
KaVo®, Dentsply Sirona®,
W&H®, Bien-Air®



- Optic
- Power 20w
- Micro head: Ø 9 x 10.8 mm

MODEL Z Micro (ultra-mini head)
CODE PA2357001
RRP £1,134

Also available in fittings for
KaVo®, Dentsply Sirona®

for regular or short shank burs

Scan/Click
For full info on
models/order codes

S-Max



OPTIC 3 + 1 FREE*

NON-OPTIC 3 FOR £1,122*

• Stainless steel • Enhanced power • 4 hole spray • DIY cartridge change • Clean head system • 2 year warranty



- Optic
- Power 26w
- Std head: Ø 12.1 x 13.3 mm

MODEL M900L
CODE P1254001
RRP £747

Also available in fittings for
KaVo®, Dentsply Sirona®,
W&H®, Bien-Air®



- Optic
- Power 23w
- Mini head: Ø 10.6x12.4mm

MODEL M800L
CODE P1253001
RRP £747

Also available in fittings
for KaVo®



- Non-Optic
- Power 26w
- Std head: Ø 12.1 x 13.3 mm

MODEL M900
CODE P1256001
RRP £696

Also available in fittings
for KaVo®, W&H®



- Non-Optic
- Power 23w
- Mini head: Ø 10.6x12.4mm

MODEL M800
CODE P1255001
RRP £696

Scan/Click
For full info on
models/order codes

Pana-Max2



1 + 1 FREE*

• Non-Optic turbine with single hole spray • Save time/money with DIY cartridge change • Clean head system • 12 month warranty



- Non-Optic
- Head: Ø 11.2 x 13.4 mm
- Coupling connection

MODEL Pana-Max2
CODE P1211001
RRP £453

for NSK PTL coupling

- Non-Optic
- Head: Ø 11.2 x 13.4 mm
- Coupling connection

MODEL Pana-Max2 KV
CODE P1212001
RRP £487

for KaVo® MULTiflex® coupling



- Non-Optic
- Head: Ø 11.2 x 13.4 mm
- Direct connection

MODEL Pana-Max2 M4
CODE P1154001
RRP £411

For Midwest M4 hose

- Non-Optic
- Head: Ø 11.2 x 13.4 mm
- Direct connection

MODEL Pana-Max2 B2
CODE P1155001
RRP £411

for Borden hose

Scan/Click
For a full list of
fittings/models/codes

*Terms and conditions apply. Offers only available through participating dealers.

NLX nano RETRO FITTED MICROMOTOR

Electric micromotors are the driving force behind speed increasing handpieces and can be fully integrated into the delivery system of most new and existing dental chairs. They can be retro fitted by an engineer or come as an option when buying new chairs. Controlling the specific speed and torque of micromotors within existing platforms is made easy thanks to built-in control pads.

Speed increasing contra-angles not only do everything a turbine does using the same FG burs, they provide a constant torque from the micromotor that does not decrease with resistance and higher load with no loss of power as they come in contact with the tooth.


Scan/Click

For more info




To discuss the offers & options contact NSK:

OFFER 1 BUY 3 Ti-Max Z95L AND GET OPTION 1 FREE



OFFER 2 BUY 4 Ti-Max Z95L AND GET OPTION 3 FREE



OFFER 3 BUY ANY OF THE BUILT-IN OPTIONS BELOW AND GET FREE Ti-Max Z95L SPEED INCREASING CONTRA-ANGLE



- option 1 NLX nano micromotor & BUILT-IN KIT *** (including MultiPad) RRP £2,229
- option 2 NLX nano micromotor & BUILT-IN KIT * + Varios170 Lux** ULTRASONIC HANDPIECE (tip options for scaling perio, endo & restorative) RRP £3,531
- option 3 NLX plus micromotor** (includes Rotary Endo functions) & BUILT-IN KIT * (including MultiPad) RRP £2,681
- option 4 NLX plus micromotor** (includes Rotary Endo functions) & BUILT-IN KIT * (including MultiPad) + **Varios170 Lux** ULTRASONIC HANDPIECE (tip options for scaling, perio, endo & restorative) RRP £3,957
- option 5 Varios170 Lux** ULTRASONIC HANDPIECE (tip options for scaling perio, endo & restorative) + BUILT-IN KIT * (including MultiPad) RRP £1,495

*Excludes install. Through participating dealers only.

NLX nano PORTABLE MICROMOTOR

- Portable electric micromotor for connecting slow, straight and speed increasing handpieces
- Puts you in control - set your own speeds
- Helps to manage aerosol generation
- Connects to your existing cart

MODEL NLX nano Portable
CODE Y1001362
RRP £2,629



OFFER NLX nano Portable ONLY £1,656 + FREE Ti-Max Z95L, WORTH £1,556



Ti-Max

2 + 1 FREE*

- Lightweight external titanium body • **Optic with internal water** • Small head/narrow neck
- **Unique spray switch on Z95L, Z45L, Z25L** • 3 year warranty



- **Optic**
 - **1:5 Increasing+Spray switch**
 - Head size $\varnothing 8.9 \times 13.4\text{mm}$
- MODEL Z95L
CODE C1034001-SW
RRP £1,556



- **Optic**
 - **1:5 Increasing**
 - 10% shorter & lighter than Z95L
 - Head size $\varnothing 8.9 \times 13.4\text{mm}$
- MODEL nano95LS
CODE C1099001
RRP £1,482



- **Optic**
 - **1:4.2 Increasing+Spray switch**
 - 45° head
 - Head size $\varnothing 9.6 \times 14.5\text{mm}$
 - Ideal for minor oral surgery too
 - For regular or extra long (25mm) FG burs
- MODEL Z45L
CODE C1064001
RRP £1,556



- **Optic**
 - **1:1 Direct Drive+Spray switch**
- MODEL Z25L
CODE C1038001-SW
RRP £1,233

Scan/Click
For a full range of
NSK Contra-angles

S-Max

3 + 1 FREE*

- Stainless steel • **All with internal water** • Choice of optic & non-optic • 2 year warranty



- **Optic**
 - **1:5 Increasing**
 - Head size $10 \times H.14.6\text{mm}$
- MODEL M95L
CODE C1023001
RRP £1,005

- **Non-Optic**
 - **1:5 Increasing**
 - Head size $10 \times H.14.6\text{mm}$
- MODEL M95
CODE C1026001
RRP £834



- **Optic**
 - **1:1 Direct Drive**
- MODEL M25L
CODE C1024001
RRP £814

- **Non-Optic**
 - **1:1 Direct Drive**
- MODEL M25
CODE C1027001
RRP £535



- **Optic**
 - **4:1 Reducing**
- MODEL M15L
CODE C1025001
RRP £952

- **Non-Optic**
 - **4:1 Reducing**
- MODEL M15
CODE C1028001
RRP £636



- **Non-Optic**
 - **1:1 Straight**
- MODEL M65
CODE H1008001
RRP £523

Scan/Click
For a full range of
NSK Contra-angles

FX

1 + 1 FREE*

- 12 month warranty



- **Non-Optic Contra-angle**
 - **1:1 Direct Drive**
 - **No internal water**
 - Push button
- MODEL FX25m
CODE C1052001
RRP £264



- **Non-Optic Straight**
 - **1:1 Direct Drive**
 - **No internal water**
- MODEL FX65m
CODE H1014001
RRP £264

NAC-EC

3 FOR £195*

- 12 month warranty



- **Non-Optic Contra-angle**
 - **1:1 Direct Drive**
 - **No internal water**
- MODEL NAC-EC
CODE Y110148
RRP £155

Scan/Click
For a full range of
NSK Contra-angles

*Terms and conditions apply. Offers only available through participating dealers.

Surgic Pro2

NEW

SURGICAL MICROMOTOR

- Compact new motor: 25% shorter, 9% lighter than current motor – Greatly improved balance
- Link by Bluetooth
- Significantly improved optics (LED) • **2 year warranty**

Complete Kit contains:

- Control unit
- Motor/Cord
- 20:1 contra angle
- Wireless foot control

- **Optic**
- **With Ti-Max X-SG20L**

MODEL Surgic Pro2
CODE Y1004195
RRP £5,856

- **Non-Optic**
- **With S-Max SG20**

MODEL Surgic Pro2
CODE Y1004196
RRP £4,164




BUY ANY Surgic Pro2 AND GET Osseo 100+ FOR FREE**



MODEL Osseo 100+
CODE Y1004176

LINK BY BLUETOOTH

- Wireless Foot Control
- iPad
- Variosurg3 (Surgical Piezo) operate with 1 wireless Foot control
- Osseo100+ confirms ISQ value over a shared interface



ADDITIONAL MOTOR/CORD



- **Optic**

MODEL SGL80M
CODE Y1004211
RRP £2,476



- **Non-Optic**

MODEL SG80M
CODE Y1004212
RRP £1,772

FREE 20:1 IMPLANT HANDPIECE OR SURGICAL STRAIGHT
WHEN YOU BUY AN ADDITIONAL MOTOR/CORD



Osseo 100+

OSSEOINTEGRATION MONITORING DEVICE

Osseo 100+ measures implant stability and osseointegration to enhance decisions about when to load the implant. This is invaluable when using protocols with shorter treatment times and treating higher risk patients. The simple steps to measure the ISQ value allow the implant loading period to be planned in advance.

MODEL Osseo 100+
CODE Y1004176
RRP £2,570

£1,799



Scan/Click
More info &
Multipeg List

*Terms and conditions apply. Offers only available through participating dealers.

Surgic Pro

SURGICAL MOTOR UNIT

Complete Kit contains:

- Control unit
- Motor/Cord
- 20:1 contra angle
- **2 year warranty**
- Corded foot control

- **Optic**
- **With Ti-Max X-SG20L**

MODEL Surgic Pro
CODE Y1003586
RRP £6,002

- **Non-Optic**
- **With S-Max SG20**

MODEL Surgic Pro
CODE Y1003587
RRP £4,481



FREE 2 ADDITIONAL 20:1 CONTRA-ANGLES*
WITH Surgic Pro COMPLETE KITS

ADDITIONAL MOTOR/CORD



- **Optic**

MODEL SGL70M
CODE E1023001
RRP £2,665



- **Non-Optic**

MODEL SG70M
CODE E1025001
RRP £1,853

IRRIGATION TUBING

- Used with Surgic Pro / Pro2 / XT and VarioSurg3



Irrigation tubing x5 (disposable)
CODE Y900079
RRP £62.20



Irrigation Tubing x1 (Autoclavable)
CODE S900055
RRP £52

Eco Friendly
autoclavable
option

FREE 20:1 IMPLANT HANDPIECE OR SURGICAL STRAIGHT

WHEN YOU BUY A SPARE MOTOR/CORD



ZYGOMATIC CONTRA-ANGLE

- **Non-Optic** • **20:1 Reduction** • **Twist chuck** • **12 Month warranty**

NSK's new handpiece for Zygomatic & Pterygoid implants. Its unique angulated shape offers optimal access and visibility.

NEW

MODEL SGX-E20R
CODE HA1200001
RRP £1,529



1 + 1 FREE* ON ZYGOMATIC CONTRA-ANGLE

IMPLANT CONTRA-ANGLES 20:1

- Max speed: 2,000min⁻¹ • **2 year warranty**



- **Optic & Dismantleable**

MODEL Ti-Max X-DSG20L CODE C1068001 RRP £1,527

- **Optic**

MODEL Ti-Max X-SG20L CODE C1003001 RRP £1,243

- **Non-Optic**

MODEL S-Max SG20 CODE C1010001 RRP £1,049

SURGICAL STRAIGHTS 1:1

- Used with standard HP burs • **2 year warranty**



- **Optic**

MODEL Ti-Max X-SG65L CODE H1009001 RRP £933

TRIPLE PACK OFFER £1,855*

- **Non-Optic**

MODEL Ti-Max X-SG65 CODE H1038001 RRP £636

TRIPLE PACK OFFER £1,308*

*Terms and conditions apply. Offers only available through participating dealers.

VarioSurg3

ULTRASONIC BONE SURGERY SYSTEM

50% OFF FOR EXISTING
Surgic Pro / Pro2 / XT USERS*†

OR

FREE 2nd HANDPIECE/CORD
OR 12 EXTRA TIPS
OR AN OSSEO 100+

- Piezo technology bone surgery system with LED.
- Lightweight slimline handpiece.
- Contains a tip stand with 6 tips – SG11, SG6D & SG7D for sinus lifts.
H-SG1 for cutting bone with SG3 & SG5 for scraping.
- Wide range of tips available for fast and more precise treatments.

MODEL VarioSurg3
CODE Y1002726
RRP £8,086

Scan/Click
For a full range of tips

Scan/Click
For the core range of tips



VarioSurg tips 1 + 1 FREE*



LINK OPTIONS

Link VarioSurg3 with Surgic Pro or Surgic Pro2

Control two units with one foot pedal. Choose a Link Stand or store everything on an iCart where you can adjust the height and move it around on its wheels.



Linking Kit for Surgic Pro

MODEL Link Stand/Cable
CODE Y1002729
RRP £512



Linking Kit for Surgic Pro2

MODEL Link Stand2
CODE Z1452001
RRP £301

MODEL Wireless Link Module
CODE Z1402001
RRP £237



iCart

MODEL iCart Duo
CODE S9090
with Surgic Pro link cable (CODE Z1189001)

MODEL iCart Duo
CODE S9090
with Surgic Pro2 link module (CODE Z1402001)


RRP £723

*Terms and conditions apply. Offers only via participating dealers. †VarioSurg3 model excludes foot pedal on 50% off promotions only, as it links to the Surgic Pro / Pro2.

VIVA ace

MOBILE DENTISTRY SYSTEM

- Ideal as a surgery back up or to visit shielded patients
- Lightweight micromotor with additional endo functions
- Multipurpose ultrasonic handpiece
- Quick setup and storage
- 4-level adjustable height and angle



MODEL VIVA ace
CODE Y1003773
RRP £12,517

Weight 8.6kg

Width 462 mm, Height 284 mm, Depth 292 mm

iProphy

NEW

CORDLESS PROPHY HANDPIECE

The complete kit comes with a 4:1 reduction head for screw-in prophy cups/brushes. Alternatively a PR-F head (optional extra) can be used that connects to disposable Prophy heads.



MODEL iProphy Complete Kit
CODE Y1500015
RRP £1,085

FREE S-Max M95L, M25L & M65 Handpieces*

OFFER 2 ×  **FREE**

Monthly **RENTAL** £433* Fully maintained for 3 years

Endodontics

Endo-Mate TC2

CORDLESS ENDO MOTOR/HANDPIECE

- Large LCD display
- 5 programmes for different file systems
- Supports most major brands of Ni-Ti files
- Torque control facility ranging from 0.1 to 4.0 Ncm depending on head:gear ratio
- Auto-reverse and alarm function sounds an audible alert when preset torque level has been reached

MODEL Endo-Mate TC2
CODE Y1001027
RRP £1,710



Endo-Mate DT2

CORDED ENDO MOTOR/HANDPIECE

- Large LCD display
- Store up to 9 speed and torque settings
- Use with Ni-Ti files from most major suppliers
- Torque control from 0.1 to 6.5 Ncm depending on head:gear ratio
- Auto-reverse on, auto-stop, autoreverse off

MODEL Endo-Mate DT2
CODE Y1004230
RRP £1,616



iPex II

APEX LOCATOR

- Accurate measurement of any root canal using SmartLogic technology
- Operates in all conditions, including dry, wet and bloody canals

MODEL iPex II
CODE Y1002208
RRP £1,015



BUY Endo-Mate DT2 AND GET iPex II FREE WORTH £1,015*

*Terms and conditions apply. Offers only available through participating dealers.

Scaler Tips


1 + 1 FREE*

TIPS FOR VARIOS
ULTRASONIC SCALER

Scan/Click
Full range
of tips

TIPS FOR AIR SCALERS

Scan/Click
Tips for
Air Scaler



AIR SCALERS



Direct connection
air scaler

- Oscillation frequency 5,600-5,850Hz
- Kit includes 3 scaler tips

MODEL AS2000
CODE T698001 (Midwest M4 hose)
CODE T697001 (Borden B2/B3 hose)
RRP £477 each



Optic air scaler with
coupling connection

- Oscillation frequency 6,200-5,850Hz
- Scratch resistant with titanium body
- Kit includes 3 scaler tips

MODEL Ti-Max S970L
CODE T1005001**
RRP £1,339

**also available in KaVo® fitting T1015001
and Sirona® fitting T1013001, RRP £1,372

FREE NAC-EC WITH AS2000

FREE 2 × NAC-EC WITH S970

PROPHY HANDPIECES

FX 1 + 1 FREE*




- Max speed 1,000 rpm
- Brushes/Cups: **Screw in**

MODEL FX57m
CODE C1054001
RRP £297

- Max speed 40,000 rpm
- Brushes/Cups: **Disposable Prophy heads**

MODEL FX65m
CODE H1014001
RRP £264



Paste defence system

NAC-EC 3 FOR £195*



AR-EC(k) AR-EC(s) NAC-EC

- Max speed 5,000 rpm
- Brushes/Cups: **Snap on cups**

MODEL AR-EC(k)
CODE Y110165
RRP £153

- Max speed 5,000 rpm
- Brushes/Cups: **Screw in**

MODEL AR-EC(s)
CODE Y110164
RRP £153

- Max speed 30,000 rpm
- Brushes/Cups: **Latch grip**

MODEL NAC-EC
CODE Y110148
RRP £155

*Terms and conditions apply. Offers only available through participating dealers.

Prophy-Mate neo SUPRAGINGIVAL AIR POLISHER



- Compact and lightweight
- Includes 60° and 80° nozzles
- 360° rotation
- Packed in a handy case
- Easily dismantled for autoclaving
- Includes 6 sachets of FLASH pearl® that quickly removes stains and plaque

CODE Y135029 for NSK PTL coupling*
 CODE Y135030 for KaVo® MULTIflex® LUX coupling
 CODE Y135031 for Sirona® quick coupling
 CODE Y135032 for W&H® Roto Quick® coupling
 CODE Y135033 for Bien-Air® Unifix® coupling

RRP £1,124



+  **FREE SECOND HANDPIECE WORTH £427***



Prophy-Mate neo nozzles

CODE T797012 60° nozzle
RRP £120

CODE T797013 80° nozzle
RRP £120

FLASHpearl

- Calcium carbonate (54µm) based powder for stain removal



4 x 300 g bottles
CODE Y900698
RRP £138

100 x 15 g sachets
CODE Y900693
RRP £189

1 + 1 FREE*

BUY ANY 2 PACKS GET A
Prophy-Mate neo NOZZLE FREE*

Perio-Mate BIOFILM ERASER



- Avoids injury to delicate periodontal pockets
- Smooth 360° rotation
- Flexible plastic nozzle tip for better access to difficult to reach areas
- Ultra-fine water spray and Perio-Mate powder eradicate biofilm in just 20 seconds
- Packed in a handy case with 160g powder

CODE Y1002654 for NSK coupling
 CODE Y1002655 for KaVo® MULTIflex® LUX coupling
 CODE Y1002656 for Sirona® quick coupling
 CODE Y1002657 for W&H® Roto Quick® coupling
 CODE Y1002658 for Bien-Air® Unifix® coupling

RRP £1,281



+  **FREE SECOND HANDPIECE WORTH £427***



- Glycine based powder (25µm)
 - 2x 160 g bottles
- Perio-Mate Powder
CODE Y900938
RRP £96.10



- Pack of 40 Perio-Mate Tip Nozzles
CODE Y1002741
RRP £96

Prophy-Mate neo + Perio-Mate

£1,854*

Varios Combi Pro iPiezo engine®

TWO ORAL CARE SYSTEMS IN ONE!



Combines ultrasonic treatments and powder therapy in one unit.
Comes with a detachable Prophy chamber. Detachable Perio chamber available as an optional extra.

MODEL Varios Combi Pro
CODE Y1002843
RRP £4,977



Portable cart
for Varios Combi Pro

MODEL iCart Prophy
CODE S9095
RRP £719

Scan/Click

FREE Perio Set WORTH £1,190*

Note: Connection to external 6mm air & water line sockets are needed to run this unit.

For more info

Varios 970 iPiezo engine®

ULTRASONIC SCALER



Ultrasonic scalers not only offer patients a more comfortable experience as less water is needed, but offer the user a wide range of 70 tips, not only for scaling/ perio/implant maintenance, but for endo, extractions, restorative & MI applications.

Units come with:


- 1 handpiece
- G4, G8 & G16 tips
- 3 tip wrenches
- sterilisation case

• **Optic**

MODEL Varios 970 LUX
CODE Y1004161
RRP £2,898

• **Non-Optic**

MODEL Varios 970
CODE Y1004160
RRP £2,494



OFFER + 2 FREE



HANDPIECES & STERILISATION CASES**

**EXCLUDING EXTRA TIPS/TIP WRENCHES

Varios2 handpiece

• **Optic**

MODEL VA2-LUX-HP
CODE E351050
RRP £531



Varios2 handpiece

• **Non-Optic**

MODEL VA2-HP
CODE E350050
RRP £279



TRIPLE PACK OFFER £1,077*

TRIPLE PACK OFFER £630*

SCALER TIPS 1 + 1 FREE*



Scan/Click
Full range
of tips

Scan/Click
To order

*Terms and conditions apply. Offers only available through participating dealers.

AFFORDABLE,
RELIABLE
& RAPID VACUUM
STERILISATION



COPPER CHAMBER

Save on running costs

- Copper chamber more efficient to run than traditional stainless steel
- Quicker instrument turnaround
- Effective drying – protects from corrosion and recontamination
- Improved reliability – reduced maintenance costs
- Extends life of equipment and handpieces



Scan/Click

For more info

**EXTENDED WARRANTY
ON NSK HANDPIECES
WITH EVERY NEW NSK
AUTOCLAVE**

iCare



Automated handpiece lubrication unit.

Compressed air connection required 5-6 bar.
6mm air line required, ideally with a shut off valve.

MODEL C2
CODE Y1002787
RRP £2,570

MODEL C3
CODE Y1002788
RRP £2,701

HALF PRICE £1,285

or **£1,065** when buying a **NSK autoclave**

RENT for just **£73.95pm[†]**

iMax S



Rapid, reliable S-type autoclave.

- 20 min fast cycles for unwrapped loads
- 25 min handpiece cycle
- Affordable, effective upgrade from non-vacuum autoclaves
- Full vacuum drying phase (Dry Handpieces!)

CODE DX91803S0200
RRP £3,224

ONLY £3,070

£350 install & delivery**

RENT for just **£115pm[†]**

iClave plus



High performance B-type vacuum autoclave.

- Easy to use
- Fast 25 min cycle for unwrapped loads
- 35 min for pouched loads
- Ideal for implant kits
- Includes Helix Test + Strips

CODE DXY1003077
RRP £6,197

FROM ONLY £3,995

£350 install & delivery**

RENT for just **£155pm[†]**

[†]Visit mysnk.co.uk/rental for more information. ^{**}For Scotland, call for details.

PROTECT IT

All-Inclusive Maintenance and Breakdown Contracts for NSK autoclaves from just **£64pm**.

One annual maintenance visit or every 1500 cycles

Additional maintenance visit for units exceeding 1500 cycles

All service parts and service sundries included

All breakdown calls*

All breakdown parts including labour*

All engineer travel in the UK*

Free on-site operator training

Free technical software upgrades

CONTACT NSK

Institutional Business Manager
Mark Beckwith
mbeckwith@nsk-uk.com 07900 246529

Ireland - Sales & Services Manager
Jonathan Singh
jsingh@nsk-uk.com 07464 675158

Scotland - Sales Manager
Angela Glasgow
glasgow@nsk-uk.com 07525 911006

North West / North East / North Wales
Clare Morris
cmorris@nsk-uk.com 07901 117669

West Midlands / South Wales
Andy Hargreaves
ahargreaves@nsk-uk.com 07872 824175

Derbys/Notts/Yorkshire/Lincs
Jo Dronfield
jdronfield@nsk-uk.com 07545 786619

Leics / Northants / CV & PE postcodes
Graham Brown
gbrown@nsk-uk.com 07541 864641

South Central / London SW, W, NW
Steve Barnes
sbarnes@nsk-uk.com 07889 181099

Central London
Teruyuki Nakanishi
tnakanishi@nsk-uk.com 07570 293317

South East & East Anglia
Francois Faro
ffaro@nsk-uk.com 07900 246822

Field Service Engineer - North
Jamie Duncombe
jduncombe@nsk-uk.com 07780 105457

Field Service Engineer - South
Lee Morris
lmorris@nsk-uk.com 07435 783707

NSK United Kingdom LTD www.mysnk.co.uk

Office 4, Gateway 1000 Arlington Business Park, Whittle Way Stevenage, SG1 2FP, UK - info@nsk-uk.com - 0800 6341909



*Offers cannot be combined. NSK reserve the right to change prices and offers without notice. Published prices are RRP's and exclude VAT and installation. Offers only available through participating dealers. Additional warranties available from NSK UK on new products and NSK UK repairs. To qualify for specific offers please send a copy invoice to NSK UK at Office 4, Gateway 1000, Arlington Business Park, Whittle Way, Stevenage, Hertfordshire SG1 2FP or email sales@nsk-uk.com. Offers end 31/12/23.

